

Austvet Endoscopy

CLEANING SOLUTIONS FOR YOUR CLINIC



\$330

OLYMPUS & PENTAX ENDOSCOPE LEAK TESTERS

Designed to detect and identify leaks from your flexible endoscopes. Stops potential fluid invasion thereby preventing the possibility of costly repairs, cross infection and downtime. **Use prior to every cleaning cycle.**

**An Economical Insurance Policy for Your Scopes!
Reusable and Easy To Use**

MEDIZYME MULTI-ENZYME SOLUTION

Medizyme Multi-Enzyme Cleaner is totally compatible with all flexible and rigid endoscopes due to its neutral PH. Provides safe and effective cleaning against proteins, carbohydrates, lipids and mucopolysaccharides.



\$15
500ml

\$110
5 Litre

PERASAFE STERILANT/DISINFECTANT

PeraSafe is a safe, fast acting disinfectant/sterilant. Supplied as a powder, PeraSafe is activated by simply dissolving in lukewarm tap water. It is a corrosion inhibitor, biodegradable and has no toxic vapour phase. Contains NO Gluteraldehyde. Suitable for use on all flexible and rigid endoscopes and accessories.



\$130
810gr Jar

ENDOSCOPE CLEANING BRUSHES


Cleaning Brush designed for maximum cleaning performance and easy handling of most instrument channels. Semi Disposable Combination Brush includes Port Stubby Brush & Channel Brush.



\$10

OLYMPUS ENDOSCOPE CLEANING DVD

Flexible Endoscope Reprocessing DVD (2008) for gastrointestinal endoscopes and bronchoscopes. This DVD provides a step-by-step demonstration of the reprocessing procedure. **Excellent reference tool.**



\$15

Austvet Endoscopy – Don't Take Risks With Your Precious Endoscopy Equipment!

www.austvetendoscopy.com.au

Prices exclude GST and Freight