

# Austvet Endoscopy



## New

### Olympus/Pentax Leakage Tester

A Great Endoscope Insurance Policy -  
Protects You, Your Patients and Your Scopes!



**\$330**

- Stop fluid potentially invading the scope
- Prevent costly repairs and downtime of your Endoscopic equipment due to moisture damage
- Decrease risk of cross infection by preventing bacteria and viruses invading the scope

*To Use: simply connect Leak Tester to ETO valve on endoscope and pressurise scope. If scope holds pressure it is water-tight, if scope loses pressure then it has a leak. Takes less than 1 minute to test!!*

*Contact Us for Further Information*

#### Other Cleaning Products for Your Endoscopes:



Medizyme Multi-Enzyme  
Cleaner

\$15 500ml  
\$110 5 Litre



PeraSafe Instrument  
Sterilant/Disinfectant

\$130  
810gr Jar



Endoscope  
Semi Disposable Brush Set

\$10



Olympus Endoscope  
Cleaning DVD

\$15