

Austvet Endoscopy

New

Olympus CX21LED Clinical Microscope

Long-Term, Cost Effective Performance for Veterinary Applications



Upright microscope with built-in LED illumination.

This imaging microscope produces a more uniform, controllable and stable illumination of samples to that provided by halogen bulbs.


\$1,500

Features:

- High-Quality clear images with a stable colour balance at any brightness level due to the incorporation of plan corrected UIS2 objectives
- LED illumination eliminates the need for frequent halogen globe replacements
- Ergonomic layout of the main controls (focus knobs, light control adjustment and stage handle)
- Mechanical stage focus lock prevents objective and specimen damage

Specifications:

- Supplied with 4x, 10x, 40x and 100x Objectives
- Expected lifetime of LED is over 20,000 hours
- Weight: 6kg



Austvet Endoscopy... Don't Take Risks With Your Precious Endoscopy Equipment!

www.austvetendoscopy.com.au

Price excludes GST and Freight