

Austvet Endoscopy

New

Oesophageal Rigid Foreign Body Removal Forceps

Extra long forceps enable foreign bodies to be manipulated back through the mouth or into the stomach.

Teeth face backwards to maximise grip and minimise damage to the oesophagus.

Simultaneous insertion of an endoscope possible in larger animals.



Working Length: 730mm
Jaw Size: 32mm L x 12mm W

\$400



Austvet Endoscopy... Don't Take Risks With Your Precious Endoscopy Equipment!

www.austvetendoscopy.com.au

Price excludes GST and Freight