

Austvet Endoscopy

QuickVet Analyser Diagnostic System



An easy to use, rapid and cost effective diagnostic unit providing fast and precise results directly in your practice.

\$3,750

Operated via a simple touch screen operation, the unit brings the “lab-on-a-chip” concept to your clinic and uses individual test cartridges based on microfluidic technology. Simply add a blood sample on the test cartridge and read the result in minutes on the 7" colour touch screen. Test for:

- **Both Activated Partial Thromboplastin Time (aPTT) and the Prothrombin Time (PT) coagulation in companion animals for detection of systemic bleeding**
- **Fibrinogen levels for horses**
- **Canine DEA1.1 Blood Typing**

Features:

- Fast results in 3-10 minutes
- Uses only 30ul of citrate blood or plasma
- Fully contained cartridge

Indications:

- Rodenticides and other toxins
- Snake bites
- Hereditary bleeding disorders
- Pre-Operative checks for at risk animals
- Systemic disease
- Disseminated Intravascular Coagulopathy (DIC)
- Blood Infusions and Typing blood donors



Cartridges Sold Separately

Austvet Endoscopy – Don't Take Risks With Your Precious Endoscopy Equipment!

www.austvetendoscopy.com.au

Prices exclude GST and Freight

