

Austvet Endoscopy


New

Askir Portable Suction Pump

A compact, high performance unit designed for both endoscopic aspiration and surgical suction



\$650



Features:

- Designed for almost continuous use up to 120 minutes
- Equipped with aspiration regulator and vacuum indicator
- Maximum suction flow 30 litres per minute
- Lightweight – only 3.6kg
- Supplied complete with autoclavable 1 litre jar with overflow valve and antibacterial filter

Austvet Endoscopy... Don't Take Risks With Your Precious Endoscopy Equipment!

www.austvetendoscopy.com.au

Price excludes GST and Freight