



## VETERINARY ENDOSCOPY SPECIALISTS

Sales • Service • Support

NEW

### Endotherapy Devices for 3 Metre Endoscopes



#### **Rat Tooth Reusable Grasping Forcep**

**\$595**

Rat Tooth Grasper offers consistent performance procedure after procedure. 15mm jaw opening. Reusable and autoclavable.

**2.4mm / 350cm Working Length**



#### **Biopsy Forcep Standard**

**\$175**

Traditional oval cup with smooth jaw. Ideal for biopsy of Equine stomach and duodenum.

**2.4mm / 350cm Working Length**



#### **Mini Oval Snare**

**\$60**

Designed with optimal stiffness to precisely capture foreign bodies. Snare opening 2.5cm

**2.3mm / 350cm Working Length**



#### **Equine Guttural Pouch Washing Catheter**

**\$30**

An endoscopically guided catheter for lung lavage or introduction into the guttural pouch to flush out the debris. Visual indicators at 2.5cm, 5cm and 7.5cm from distal tip to aid in introduction. Sterile.

**2.3mm Diameter / 350cm Working Length**



#### **Semi Disposable Channel Cleaning Brush**

**\$15**

For cleaning flexible endoscope channels between 2mm and 4.2mm. **350cm Working Length.**

Austvet Endoscopy - Don't Take Risks With Your Precious Endoscopy Equipment!

[www.austvetendoscopy.com.au](http://www.austvetendoscopy.com.au)

\*Prices exclude GST and Freight\*