



AUSTVET ENDOSCOPY

PRODUCT CATALOGUE 2023

Est. 2004

OLYMPUS®

True to Life



Make the Invisible, Visible!

Don't Take Risks With Your Precious Endoscopy Equipment!

Please don't hesitate to call any time with any questions or queries
or come and visit our showroom soon:

Unit 5 / 14 Lionel Rd, Mt Waverley VIC 3149

Maria: 0409 383 327 | Ennio: 0408 053 100 | Wally: 0438 084 418

sales@austvetendoscopy.com.au www.austvetendoscopy.com.au www.steriflex.com.au



| | |
|------------------------------------|----|
| Olympus Partnership | 3 |
| Olympus HD Video Endoscopy System | 4 |
| Endoscopic Auxiliaries | 5 |
| Olympus Flexible Video Endoscopes | 6 |
| Portable HD Video Endoscopy System | 8 |
| Portable Endoscope | 9 |
| Retrieval Devices | 10 |
| Therapeutic Devices | 11 |
| Fluid Administration & Therapy | 12 |
| Endoscopic Accessories | 13 |
| Arthroscopic Fluid Management | 15 |
| Otoscopy/Rhinology | 16 |
| Portable HD Video System | 17 |
| Arthroscopy | 18 |
| Laparoscopy | 24 |
| Caiman Vessel Sealing Device | 27 |
| Trocars & Cannulas | 28 |
| EO Gas Sterilisation | 32 |
| Service & Repair | 34 |
| Cleaning Solutions | 35 |

DISCLAIMER: Please note that prices and availability are correct at the time of publishing but may be subject to change without notice.

***All prices exclude GST and freight**


OLYMPUS®
True to Life

EVIS X1 Video System Centre CV-1500
Elevating the Standard of Endoscopy

Olympus Australia's corporate philosophy is to make people and their much loved animal's lives healthier, safer and more fulfilled. One of their purposes is to *make the invisible, visible to see and treat disease*. They are on the cutting-edge of medical procedures and scientific discoveries and deliver a range of quality services and State of the Art endoscopy products.

The philosophy, "True to Life" sits in the heart of their beliefs and every day they develop and deliver medical technological solutions that elevate the standard of care to see and treat diseases with a full conviction that every life counts.

The medical instrument division encompasses a broad range of endoscopy and surgical instruments, with particular emphasis in the fields of:

Gastroscopy, colonoscopy, bronchoscopy and the entire spectrum of endoscopic and keyhole surgery, including laparoscopy, cystoscopy and arthroscopy.

Austvet Endoscopy is the exclusive veterinary distributor of Olympus Endoscopy products in Australia and New Zealand.

Recognised as the world's leading supplier of medical endoscopy equipment, we are proud to offer the Olympus range of products to the veterinary industry and support their values and philosophies.



OLYMPUS®

The Olympus 190 Video System is cutting-edge technology for unsurpassable visualisation.

Endoscopy Video System

Includes:

- Olympus CV-190 Video Processor with Keyboard and Scope Cable (MAJ-1430)
- Olympus CLV-190 Xenon Light Source with Water Bottle
- Endoscopy Trolley with Olympus OEV-261H High Definition 26" Monitor, Fully Adjustable Monitor Arm and Endoscope Hanger

Features:

- Full HD imaging capabilities when used in conjunction with H-Series video endoscopes.
- Narrow Band Imaging to enhance the visibility and characterisation of healthy and abnormal tissue and other structures on the mucosal surface
- Compatible with all Olympus Endoscopes & Surgical Camera Heads
- Image capture system is optional (see page 5)
- Video Camera head is optional

\$29,950



Camera Head Optional Extra



**New Endoscopy Trolley****\$1,950**

Trolley designed especially for Video Endoscopy.

Features:

- Heavy duty steel construction
- 4 tiers
- Medical grade castors with 100mm wheels
- Front castors with brakes
- Conceals all electrical wiring
- Power board 6-way, fully adjustable monitor arm and endoscope hanger

Dimensions: 550mm W x 550mm D x 1350mm H
(Height excludes monitor)

Monitor Optional Extra**New Microsoft Surface Laptop 4****\$3,490**Microsoft Surface Laptop 4 with custom **HD** Image Capture Software.**Features:**

- 13.5" display
- 2256 x 1504 screen resolution
- **HD** image and **HD** video capture capabilities
- Ability to input patient details (up to 80 characters)
- Supplied with SDI or DVI to USB converter

Also note: Due to computer technology changing rapidly specifications and price may differ. Please enquire.



Specifications & Compatibility Chart

| Endoscope Model | Working Length | Diameter | Biopsy Channel | Angulation |
|-----------------|----------------|----------|----------------|------------|
| URF-V | 670mm | 3.2mm | 1.2mm | 2 Way |
| URF-V2 | 670mm | 2.8mm | 1.2mm | 2 Way |
| CYF-V2 | 380mm | 5.4mm | 2.2mm | 2 Way |
| CYF-VH | 380mm | 5.5mm | 2.2mm | 2 Way |
| BF-P180 | 600mm | 4.9mm | 2mm | 2 Way |
| BF-Q180 | 600mm | 5.2mm | 2mm | 2 Way |
| BF-P190 | 600mm | 4.2mm | 2mm | 2 Way |
| BF-Q190 | 600mm | 4.8mm | 2mm | 2 Way |
| BF-H190 | 600mm | 5.5mm | 2mm | 2 Way |
| GIF-XP150N | 1100mm | 5.5mm | 2mm | 4 Way |
| GIF-N180 | 1100mm | 4.9mm | 2mm | 2 Way |
| GIF-XP180N | 1100mm | 5.5mm | 2mm | 4 Way |
| GIF-Q180 | 1030mm | 8.8mm | 2.8mm | 4 Way |
| GIF-H180 | 1030mm | 9.8mm | 2.8mm | 4 Way |
| GIF-H190 | 1030mm | 9.2mm | 2.8mm | 4 Way |
| PCF-H180AL | 1680mm | 11.8mm | 3.2mm | 4 Way |
| PCF-H190AL | 1680mm | 11.7mm | 3.2mm | 4 Way |

**All Second Hand Endoscopes Pending Availability*



| Image | Air / Water Facility | Processor CV-190 | Processor CV-180 | Price Second Hand |
|-------|----------------------|------------------|------------------|-------------------|
| | No | ✓ | ✓ | \$8,900 |
| | No | ✓ | ✗ | \$11,900 |
| | No | ✓ | ✓ | \$3,500 |
| | No | ✓ | ✓ | \$6,900 |
| | No | ✓ | ✓ | \$8,900 |
| | No | ✓ | ✓ | \$8,900 |
| | No | ✓ | ✗ | \$12,900 |
| | No | ✓ | ✗ | \$12,900 |
| | No | ✓ | ✗ | \$13,500 |
| | Yes | ✓ | ✓ | \$9,500 |
| | Yes | ✓ | ✓ | \$11,900 |
| | Yes | ✓ | ✓ | \$13,900 |
| | Yes | ✓ | ✓ | \$7,900 |
| | Yes | ✓ | ✓ | \$10,900 |
| | Yes | ✓ | ✗ | \$14,900 |
| | Yes | ✓ | ✓ | \$10,900 |
| | Yes | ✓ | ✗ | \$14,900 |



- Dimensions: 48cm x 42cm x 29cm
- Weight: 17kg

Portable battery powered **HD** video-endoscopy system designed for ease of use for Veterinary surgeons performing upper airway and gastroscopic examinations.

Endoscopy Video System

\$16,500

Includes:

- Battery Powered Video Processor with battery life indicator (over 3 hours of battery life)
- Microsoft Surface Book with 13.5" Screen
- Custom Made Carry Case
- Easy to use **HD** Image/Video Capture Software

VET-9215 HD Video Endoscope

\$15,900

Specifications:

- Working Length: 1500mm
- Diameter: 9.2mm
- Biopsy Channel: 2.8mm
- Air/Water Facility
- 4 Way Angulation
- Includes Leak Tester, Water Bottle and Cleaning Accessories

VET-6010 HD Video Endoscope

\$16,500

Specifications:

- Working Length: 1000mm
- Diameter: 6mm
- Biopsy Channel: 2mm
- Air/Water Facility
- 4 Way Angulation
- Includes Leak Tester, Water Bottle and Cleaning Accessories

VET-9830 HD Video Endoscope

\$16,500

Specifications:

- Working Length: 3000mm
- Diameter: 9.8mm
- Biopsy Channel: 2.8mm
- Air/Water Facility
- 4 Way Angulation
- Includes Leak Tester, Water Bottle and Cleaning Accessories

VET-3512 Video Endoscope

\$15,690

Specifications:

- Working Length: 600mm
- Diameter: 3.5mm
- Biopsy Channel: 1.2mm
- 2 Way Angulation
- Includes Leak Tester, Water Bottle and Cleaning Accessories





AOHUA MVE-9215 Portable Endoscope

When portability is A MUST

- Working Length: 150cm
- Diameter: 9.2mm
- Biopsy Channel: 2.8mm
- 4 Way Angulation
- LED Illumination

Scope with Small Monitor

\$11,250



Clax Smart Fold Cart

\$490

- Can be collapsed quickly and easily at the push of a button
- Has two large platforms, supplied standard with one collapsible storage box
- Has many useful features such as foot brake, holding device for boxes, upward folding platform and removable wheels for space-saving storage
- Made of injection moulded plastic and aluminium for long lasting use




Specifications:

| | |
|--|--|
| 60kg Total capacity | 7kg Unit weight |
| 20kg Top load capacity | 40kg Bottom load capacity |
| Folding Handle type | 1 crate Included |
| 890 x 550 x 1030mm Dimensions assembled | 670 x 470 x 110mm Dimensions folded |

**Great for transporting the Portable
Equine Endoscopy System**

New Oesophageal Rigid Foreign Body Removal Forceps

Teeth face backwards to maximise grip and minimise damage to the oesophagus.

| | Diameter at Jaws in mm | Length at Jaws in mm | Shaft Diameter in mm | Shaft Length in mm | Price |
|---|------------------------|----------------------|----------------------|--------------------|--------------|
|  | 13 x 9 | 34 | 7.5 x 6 | 640 | \$500 |
|  | 11.5 x 6 | 31 | 6.5 x 5 | 500 | \$500 |
|  | 5.6 x 4.2 | 15.5 | 4 x 3 | 400 | \$500 |



"Standard" Size Most Popular

New Flexible Retrieval Devices**Retrieval Grasping Basket**

For the extraction of round, slippery objects.

Paediatric: For 2.0mm > Instrument Channels, 180cm Length Basket Opening 32mm

Adult: For 2.8mm > Instrument Channels, 230cm Length Basket Opening 30mm

\$300

Special order for finer instrument channels < 1.2mm

\$470

Retrieval Tripod Prongs

For the collection of irregular shaped objects.

Paediatric: For 2.0mm > Instrument Channels, 180cm Length

Adult: For 2.8mm > Instrument Channels, 230cm Length

\$150

SnareMaster Oval Loop

Designed with optimal stiffness to precisely capture foreign bodies. Handle shaft is marked at 5mm increments to indicate the snare throw opening.

Adult: For 2.8mm > Instrument Channels, 230cm Length Snare Opening 30mm

For 2.8mm > Instrument Channels, 230cm Length Snare Opening 125mm (Jumbo)

For 2.8mm > Instrument Channels, 350cm Length Snare Opening 25mm

\$50

SnareMaster Crescent Shape

Designed with optimal stiffness to precisely capture foreign bodies. Handle shaft is marked at 5mm increments to indicate the snare throw opening.

Paediatric: For 2.0mm > Instrument Channels, 230cm Length Snare Opening 25mm

\$50

A Cord for Hot Snare

\$80

Spider-Net Retrieval Device

For the collection of round, blunt or otherwise hard to grasp foreign bodies. The "spider" net delivers a flexible, durable pouch eliminating the risk of accidental loss during retrieval.

Adult: For 2.8mm > Instrument Channels, 230cm Length Net Opening 3x6cm

\$150

Ratigator Reusable Grasping Forceps

Combines a rat tooth with alligator jaws into one instrument. 15mm jaw opening. Reusable & autoclavable.

Paediatric: For 2.0mm > Instrument Channels, 210cm Length

Adult: For 2.8mm > Instrument Channels, 230cm Length

For 2.8mm > Instrument Channels, 350cm Length

\$650

Foreign Body Removal Hood

\$80

This hood features a soft bell shaped hood for the removal of sharp or pointed gastric foreign bodies and folds backwards on insertion. The hood prevents exposure of the oesophageal and pharyngeal wall to injuries from the edges of the foreign body. Latex Free.





New Flexible Biopsy Forceps

| | Instrument Channel | Working Length | Price |
|--|--------------------|----------------|--------------|
| Biopsy Forceps/Rat Tooth (reusable) | 1.2mm | 115cm | \$950 |
| Serrated Cup Paediatric | 2.0mm | 160cm | \$50 |
| Oval Cup Paediatric | 2.0mm | 160cm | \$20 |
| Serrated Cup w Spike Paediatric | 2.0mm | 160cm | \$50 |
| Serrated Cup Adult | 2.8mm | 230cm | \$20 |
| Serrated Cup w Spike Adult | 2.8mm | 230cm | \$20 |
| Oval Cup w Spike Adult | 2.8mm | 230cm | \$20 |
| Oval Cup Adult | 2.8mm | 230cm | \$20 |
| Oval Cup Adult | 2.8mm | 400cm | \$100 |



Flexible Cytology Brushes

Designed to enhance the efficiency of bronchoscopic and gastrointestinal examinations and improve collection of cytology samples. Features optimal brush hardness for better sample collection and thinner bristles for improved cell collection.

| | Instrument Channel | Working Length | Brush Size | Price |
|-------------------------|--------------------|----------------|------------|-------------|
| Bronchial | 1.2mm | 115cm | 2mm | \$50 |
| | 2.0mm | 120cm | 3mm | \$30 |
| | 2.0mm | 120cm | 2mm | \$30 |
| Gastrointestinal | 2.0mm | 200cm | 3mm | \$30 |
| | 2.8mm | 200cm | 3mm | \$30 |



Injector Needle

Designed to provide proper flexibility and rigidity to ensure safe piercing and excellent puncture performance. Ideal for injecting Formalin to treat Ethmoid Haematomas in horses. Assorted needle lengths and gauges. Large range of needle gauges and lengths for paediatric and adult scopes, please enquire

| Instrument Channel | Working Length | Needle Diameter | Needle Length | Price |
|--------------------|----------------|-----------------|---------------|-------------|
| 2.0mm | 200cm | 23G | 5mm | \$70 |
| 2.8mm | 230cm | 23G | 6mm | \$70 |
| 2.8mm | 240cm | 21G | 8mm | \$70 |



Olympus EZDilate Inflation Device

\$150

Used while performing balloon dilation procedures to inflate the balloon, monitor the pressure within the balloon and deflate the balloon. Easy to read luminescent dial and quick release lock lever mechanism. Maximum pressure of 30 atmospheres.




Olympus EZDilate Multi-Stage Endoscopic Balloon Dilators

\$350

Designed for oesophageal dilation under direct endoscopic vision. Low compliance balloon material produces dependable, controllable radial dilation. Low profile balloon design assures smooth accessibility into difficult strictures.

Assorted Sizes: For 2.8mm > Instrument Channels, from 6mm to 20mm Diameter

| Inflation diameter (mm) | 6, 7, 8 | 8.5, 9.5, 10.5 | 11, 12, 13 | 13.5, 14.5, 15.5 | 16, 17, 18 | 18, 19, 20 |
|--------------------------|-------------|----------------|-------------|------------------|-------------|-------------|
| Inflation pressure (atm) | 2, 5.5, 9.5 | 2, 5.5, 9.5 | 2, 3.5, 5.5 | 2.5, 4, 6 | 2.5, 4, 5.5 | 2.5, 3.5, 5 |
| Balloon length (mm) | 65 | 65 | 55 | 55 | 55 | 55 |



Mila Double/Triple Lumen Flushing Catheter

Flushing catheter for guttural pouch lavage. Can also be used as a sampling catheter for sterile specimen retrieval via an endoscope. The set contains a glycol plug in the outer catheter and an inner catheter for retrieval of the sterile specimen.

1.8mm Diameter/250cm Length/Double Lumen
2.3mm Diameter/190cm Length/Triple Lumen

\$80
\$80



Equine Uterine Flushing Catheter

Used for flushing and irrigation, administration of medications, embryo collection and transfer. Autoclavable. Sterile.

12mm Diameter/80cm or 150cm Length

\$150



Endoscopic BAL Small Animal/Equine Single Lumen Flushing Catheter

An endoscopically guided catheter for lung lavage or introduction into the guttural pouch to flush out pouch debris. Visual indicators at 2.5cm, 5cm and 7.5cm from the distal tip to aid in introduction. Sterile.

1.8mm Diameter/120cm Length
1.8mm Diameter/250cm Length
2.3mm Diameter/220cm Length
2.3mm Diameter/400cm Length

\$30
\$50
\$15
\$20



Silicone Stomach Tube

Three strategically placed fenestrations at the distal tip provide improved safety and effectiveness. Includes syringe and funnel adaptor. Autoclavable. Non-Sterile. (Internal dimensions listed).

4.5mm/91cm Length or 5.0mm/91cm Length
5.0mm/183cm Length or 6.5mm/150cm Length
6.5mm/300cm Length or 9.5mm/300cm Length
11mm/300cm Length or 16mm/305cm Length
10mm/300cm Length (Standard Line)

\$80
\$80
\$115
\$150
\$70



Broncho-Alveolar Lavage Catheter (Silicone) BAL TUBE

The BAL catheter is used for lung lavage of selected segments of the lung fields in horses. Equipped with a distal balloon that can be inflated to prevent retrograde fluid flow and a proximal two-way fitting allows simultaneous connection for both irrigation and aspiration. Sterile and Autoclavable

(The plastic proximal two-way fitting and three-way stopcock are not autoclavable).

10mm Diameter/240cm Length
10mm Diameter/300cm Length
Three-Way Replacement Stopcock

\$170
\$180
\$20





Standard Single Scope Hanger

Requires only two screws for mounting.
Dimensions: 145mm W x 135mm D

\$150



Olympus, Aohua and Pentax Water Bottles

Pentax Water Bottle to suit all Pentax Endoscopes
Olympus Water Bottle Autoclavable for 140/160/180 Series
Olympus Endoscopes
Olympus Water Bottle Autoclavable for 100 Series
Olympus Endoscopes

\$500

\$280

\$190

Second Hand Bottles Also Sold Pending Availability



Olympus, Aohua and/or Pentax Leakage Testers

Designed to detect and identify leaks from your flexible endoscopes. Stops potential fluid invasion thereby preventing the possibility of costly repairs, cross infection and downtime.

\$300

ALL flexible endoscopes MUST be leak tested prior to EVERY submersion to ensure YOU, YOUR PATIENTS AND YOUR SCOPES are protected.



Askir Endoscopic Suction Pump & Accessories

Compact, high performance unit designed for both endoscopic aspiration and surgical suction. Unit is supplied complete with polycarbonate autoclavable 1 litre jar with overflow valve and antibacterial filter. Pump is equipped with aspiration regulator and vacuum indicator. Designed for almost continuous use up to 120 minutes. Maximum suction flow 30 litres per minute. Portable and lightweight.

\$690

Accessories:



Mucus Specimen Trap
option of 40 cc and 70cc

\$10



Suction Handpiece Frazier
2.6mm, 3.3mm or 4.0mm

\$25



Suction Handpiece Yankauer 29cm
Standard or Paediatric

\$35





NEW AE-L100 Endoscopic LED Light Source (No Air Pump)

\$1,350

Ideal for use with all rigid endoscopes including Arthroscopes, Laparoscopes, Oscopes, Rhinoscopes and Cystoscopes.

Features:

- Extremely compact (30cm x 20cm x 10cm)
- Expected lifetime of LED is over 40,000 hours
- Extremely bright light
- Supplied with light guide cable



For Rigid Endoscopes

LED Portable High-Intensity Light Handles for Rigid and Flexible Endoscopes Now with Variable Brightness

LED Portable Light Handles are an easy option for quick examinations in the consult room or for fieldwork. They feature high-intensity illumination and are a great light source backup. Handles attach to most rigid and flexible endoscopes.



For Flexible Endoscopes

Features:

- Lithium-Ion rechargeable battery
- Extremely compact
- High-intensity white light with variable brightness

Kit Includes:

- LED light source handle
- Rechargeable Lithium-Ion battery
- Battery charger

LED Light Handle for Rigid Endoscopes **\$250**

LED Light Handle for Flexible Endoscopes **\$300**



New Sterile Camera Sleeve with Telescope Adaptor

Disposable surgical sterile video camera sleeve/cover with enclosed adaptor for telescope connection. Safe, quick and easy application. Latex free. Offers 100% sterility without chemicals to reduce risk of infection. Excellent equipment protection. Dimensions: 14cm W x 285cm L

\$30 Each
or
\$1,080
Box of 40



Olympus UHI-3 35 Litre Insufflator

\$690

Olympus insufflator which features preset pressure and flow rate, leakage warning display and overpressure warning display and alarm. Supplied with hose for gas bottle. Designed to support a wide range of Laparoscopic insufflation applications.



New Laparoscopic Insufflator Tubing

\$10

Disposable CO2 Insufflator Filter with 3m Tubing. Has universal fitting however any additional connectors required can be supplied at no charge. Filtration area 24cm² for added protection and greater gas flow. Sterile and Latex free.



New Verres Needles

Verres Needles equipped with a spring loaded obturator to insufflate the abdomen in Laparoscopic surgery.

| | |
|-------|--------------|
| 100mm | \$145 |
| 150mm | \$163 |
| 200mm | \$210 |



New Linvatec 10k Pump. The Simplest High-Performance Pump in Arthroscopy.

Ease of Use

- Intuitive one-step/one-hand cassette loading
- Quick Setup and priming
- Integrated pole mount works with any standard I.V. pole
- Compatible with video cart shelf or any flat surface

Superior performance

- Immediate response
- High flow
- Maintains distension under even the most demanding situations
- Robust autoclavable remote control

Specifications

- Maximum set pressure 150mmHg
- Maximum flow rating 2000ml/min
- Weight: 3kg
- Dimensions: 16.5cm tall, 20.3cm deep and 19.7cm wide
- Voltage 90–240v

\$5,500



10k100 Endoscopy Tube Set

\$700

box of 10



New Dr Fritz Pressure Infusion Cuff

| | | |
|----------------------------------|-------------------|--------------|
| 500ml - 1 Litre (with manometer) | AR-08-100N | \$450 |
| 3 - 5 Litre (with manometer) | AR-08-300 | \$750 |



**New Standard Veterinary Rhinoscope Set****\$5,750**

Rhinoscope set supplied with two sheaths, with the larger sheath allowing both biopsy and irrigation.

Features universal eyepiece for compatibility with all Endoscopic video cameras.

Multi-purpose rigid endoscope which can also be used for:

Arthroscopy

Avian Endoscopy

Cystoscopy

Laparoscopy

Thoracoscopy

Vaginoscopy

System Includes:

- Rhinoscope 2.7mm 30 Deg (190mm Working Length)
- 3.7mm/11 Fr Sheath (no biopsy facility – only irrigation)
- 4.7mm/14 Fr Sheath (biopsy & irrigation)
- High-intensity LED Portable Light Handle (with variable brightness)
- Cleaning Brush & Custom Carry Case
- Supplied with your choice of 1.6mm Biopsy or Grasping Forceps (see page 17)

Otoscope not included (please refer to below)

New Standard Veterinary Otoscope Set**\$1,990**

Complete, all-in-one veterinary otoscope. Ideal for quick examinations in the consult room.

Large field of view and depth of field.

Also used for foreign body removal and/or biopsy and quick diagnostics.

Set Includes:

- Otoscope 5.5mm 0 Deg (100mm Working Length)
- High-intensity LED Portable Light Handle (with variable brightness)
- Two Way Stopcock (Removable)
- Cleaning Brush & Carry Case (refer to image of carry case on page 16)
- Supplied with your choice of 1.6mm Biopsy or Grasping Forceps (see page 17)



Essential Instrument for Every Clinic



New All-in-One Portable Video System (AE-1080 HD)

\$7,950

This all-in-one portable **HD** Video System is a complete, user friendly solution for all types of rigid endoscopy for the practice as well as the mobile vet.

Features:

- 15.5" full **HD** display
- 3 metre long camera and light lead cables
- Image/video capture facilitated from camera head
- Multiple video outputs for slave monitors

System Includes:

- 15.5" full **HD** display
- Video processor
- Camera head
- LED light source
- Light Lead
- Custom carry case

New 1.6mm/5fr Semi Flexible Biopsy/Retrieval Forceps

\$350

Will suit rigid ENT/Otoscopes with instrument channels 1.6mm and above.

Working length: 270mm.



48404

Ultra slim 1.6mm/5fr semi flexible biopsy forcep. Ideal for biopsy of ear, nose or throat.



48405

Ultra slim 1.6mm/5fr semi flexible grasping forcep. Ideal for removal of foreign bodies from ear canals and nasal cavities.



48405-B

Ultra slim 1.6mm/5fr semi flexible dual purpose biopsy/retrieval forcep.

New Olympus A37026A ENT/Cystoscopy Endoscope & Accessories

\$9,520

Compact slim design and lightweight, this scope is the ideal choice for ENT and Cystoscopy. Integrated stopcock saves time and simplifies procedure. Atraumatic smooth beak ensures optimum patient comfort.

Supplied with instrument tray and sealing caps.

Specifications:

- 7 Degree field of View
- 2.6mm Diameter
- 160mm Working Length
- 1.4mm Instrument Channel
- 45 Degree Angled Ocular
- Fully Autoclavable

Accessories:

New Flexible Grasping Forceps
1mm / 280mm Working Length

\$700


New Flexible Biopsy Forceps
1mm / 280mm Working Length

\$600



2.7mm

A70963A Telescope,
\$3,090 Autoclavable,
 30° direction of view,
 Wide angle field of view,
 70 mm working length



WA70992A Trocar tube,
\$1,090 Rotatable,
 1 Stopcock,
 Standard flow,
 for 70 mm telescopes



A70978A Trocar spike,
\$270 triangular,
 for 70 mm telescopes



A70979A Trocar spike,
\$270 blunt,
 for 70 mm telescopes



A00507A Sealing, set
\$70 Set of spare sealing rings:
 Contains 2 large sealings
 (Transparent) and 5 small sealings (black)



A0460 Light-guide adapter,
\$70 for OES 4000 light-guide cables

00332 Light-guide adapter,
\$25 for Wolf light-guide cables


A5333 Eyepiece cup,
\$105 for telescopes





A70941A
\$3,400

Telescope,
Autoclavable,
30° direction of view,
115° field of view,
160mm working length



WA70952A
\$1,090

Trocar tube,
Rotatable,
1 Stopcock,
Standard flow



A70954A
\$270

Trocar spike,
Triangular




WA70953A
\$1,090

Trocar tube,
Rotatable,
2 Stopcocks,
Standard flow



A70955A
\$270

Trocar spike,
Blunt




New Standard Line 3mm ENT/Arthroscope

Arthroscope supplied with 4mm Cannula and Sharp and Blunt Trocar.

\$1,650

Specifications:

- Working Length: 175mm (with cannula 135mm)
- Diameter: 3mm
- Direction of View: 30 Degree



New Standard Line 4mm ENT/Arthroscope

Arthroscope supplied with 5mm Cannula and Sharp and Blunt Trocar.

\$1,650

Specifications:

- Working Length: 175mm (with cannula 135mm)
- Diameter: 4mm
- Direction of View: 30 Degree



New Dr Fritz Small Animal Arthroscopy Hand Instruments

Instruments provide easy handling and an ergonomic design with handles. All instruments are specifically designed for Small Animal Arthroscopy.



AR-17-03S

Bayonet Knife with Handle, 2.7mm

\$390



AR-17-05S

Banana Knife with Handle, 2.7mm

\$390



AR-17-12S

Mikro Picking Knife with Handle, 2.7mm

\$390



AR-17-19S

Hook Knife with Handle, 2.7mm

\$390



AR-17-06S

Curette Curved with Handle, 2.7mm

\$390



AR-17-07S

Ring Curette with Handle, 2.7mm

\$390



AR-17-11S

Meniscutome Smilie with Handle, 2.7mm

\$390



AR-17-09S

Hook Probe with Handle, 2.2mm

\$390



AR-17-13

Special Elevator with Handle, 2.7mm (Metal)

\$390



AR-17-22

Coronoid Chisel with Handle, 3/4mm Working Length 20cm (Metal)

\$390



I-17-822 (SET)

Milling Drill 2.2mm, Working Length 12cm, Rotating Handpiece and Long Lasting Drill Steel (Dismantable)

\$550



Replacement Drill Steel (for Milling Drill)

\$180



S110-1935

Protection & Transportation Sheath for Small Arthroscopes 12cm, Autoclavable

\$45

New Dr Fritz Small Animal Arthroscopy Punches, Forceps & Scissors

AS-17-521

\$1,550

Mini Scissors for Arthroscopy, 2mm,
Working Length 11cm

AT-19-231

\$1,550

Mini Rongeur & Biopsy Forceps, 2.2/2.7mm,
Working Length 11cm

AS-17-632

\$1,550

Universal Rongeur & Grasping Forceps, 2.7mm,
Working Length 12cm

AS-17-662

\$1,550

Alligator Grasping Forceps with Ratchet, 2.7mm,
Working Length 12cm

AS-17-621

\$1,550

Arthroscopic Hook Scissor, 2.7mm,
Working Length 12cm

AS-17-672

\$1,550

Arthroscopic Punch, 2.7mm,
Working Length 12cm

AS-17-641

\$1,550

Biopsy & Spoon Forceps, 2.7mm,
Working Length 12cm

New Dr Fritz Small Animal Working Cannula Set

This system allows the specific exchange of different working cannulas. Beginning with the smallest diameter, the vet can always change to a bigger diameter by using the changing rods supplied. Cannulas are specifically designed for Small Animal Arthroscopy.

I-17-755

\$990

Working Cannula Set

Consists of:

- 3 Cannulas for Elbow 2.3mm, 2.9mm and 3.5mm, Working Length 3cm – 5.5cm
- 2 Cannulas for Shoulder 2.9mm and 3.5mm, Working Length 7cm – 8cm
- 2 Sharp Trocars
- 2 Changing Rods, 2.3mm and 2.8mm



New Dr Fritz Equine Arthroscopy Hand Instruments

Instruments have been specifically designed for Equine Arthroscopy.

| | | |
|---|--|--------------|
|  | I-18-100 Osteotome Curved with Handle, 4mm | \$240 |
|  | I-18-110 Osteotome Straight with Handle, 4mm | \$240 |
|  | I-18-520 Hook Probe with Handle, 3mm, Working Length 12cm | \$240 |
|  | I-18-521 Special Knife with Handle, 3mm, Break Safe Special Cut | \$285 |
|  | I-18-963 Hook Scissors, Straight, Single Action Jaws, 4.5mm, Working Length 18cm | \$840 |
|  | I-18-601 Curette Straight with Handle, Oval Cup 6 x 4mm | \$240 |
|  | I-18-621 Curette, 35 Degree Oblique Left, Oval Cup, 6 x 4mm | \$240 |
|  | I-18-622 Curette, 35 Degree Oblique Right, Oval Cup, 6 x 4mm | \$240 |
|  | I-18-518 Egress Cannula with Sharp Trocar and Stopcock, 2.7mm, Working Length 7cm | \$295 |
|  | I-18-519 Egress Cannula with Sharp Trocar and Stopcock, 4mm, Working Length 7cm | \$295 |
|  | I-18-635 Special Egress Curette Straight, Oval Cup, 5.5 x 3mm and Leur Lock (f) Suction Connector, Working Length 12cm | \$450 |
|  | I-18-536 Special Elevator for Coffin Joint, Working Length 8cm | \$260 |
|  | I-18-517 Elevator Curved with Handle, 6mm | \$300 |
|  | I-18-534 Special Elevator for Sesamoide Fractures, Curved, Long 4mm | \$350 |



New Dr Fritz Equine Arthroscopy Line Instruments

“Standard Line” Rongeur

I-18-6611

\$1,050

Straight, Fenestrated, Oval Cup 4 x 8mm,
Working Length 11cm

I-18-662

\$1,050

Angled Upwards, Fenestrated, Oval Cup 4 x 8mm,
Working Length 11cm



“Power Line” Rongeur

I-18-973

\$1,290

Straight, Heavy Durable Handle, Oval Cup 5 x 10mm,
Working Length 15cm

I-18-974

\$1,290

Straight, Heavy Durable Handle, Fenestrated, Small Tip,
Oval Cup 5 x 10mm, Working Length 15cm



New Dr Fritz Equine Arthroscopy Ergonomic Instruments

Straight Rongeur

I-18-7613

\$1,290

Straight, Ergonomic Handle, Fenestrated
Cup 4 x 8mm, Working Length 12cm



Hook Grasper “Chip Catcher” Rongeur

I-18-6414

\$1,290

For Large Fragments, Ratchet,
Working Length 12cm



New Ferris Smith Equine Rongeurs

Premium quality German made Rongeurs.
Features straight and angled upward tip
with 180mm working length.
Lifetime warranty.

\$1,090

Sizes Available:

- 3mm W x 10mm L Tip (Straight & Angled Upward)
- 4mm W x 10mm L Tip (Straight & Angled Upward)
- 5mm W x 10mm L Tip (Straight & Angled Upward)
- 6mm W x 10mm L Tip (Straight only)



New Dr Fritz Needle Holder with Handle

Needle self positioning, axial handle.

Small Animal: **L-19-562V**

Left Turned 5mm,
Working Length 22cm

\$1,150

L-19-563V

Angled Upwards 5mm,
Working Length 22cm

\$1,150

New Dr Fritz Small Animal 5mm Hand Instruments

Instruments are easily dismantable and rotatable through 360 Degrees. Each Instrument is supplied complete with jaw insert, sheath and ratchet handle.



531.61H01V

Biopsy Forceps, Double Action Jaws, With Spikes
5mm, Working Length 33cm

\$825

541.61H01V

Biopsy Forceps, Double Action Jaws, With Spikes
5mm, Working Length 43cm

\$825



531.01H05V

Mini Metzenbaum Scissors, Curved Left, Insulated
5mm, Working Length 33cm

\$825

541.01H05V

Mini Metzenbaum Scissors, Curved Left, Insulated
5mm, Working Length 43cm

\$825



531.22H11V

Grasping Forceps Universal, Insulated
5mm, Working Length 33cm

\$825

541.22H11V

Grasping Forceps Universal, Insulated
5mm, Working Length 43cm

\$825



531.40H01V

Dissecting Forceps Maryland, Turned, Insulated
5mm, Working Length 33cm

\$825

541.40H01V

Dissecting Forceps Maryland, Turned, Insulated
5mm, Working Length 43cm

\$825



531.30H12V

Grasping Forceps Claw, Straight, Metal
5mm, Working Length 33cm

\$825

541.30H12V

Grasping Forceps Claw, Straight, Metal
5mm, Working Length 43cm

\$825



531.27H11V

Grasping Forceps Babcock
5mm, Working Length 33cm

\$825



New Dr Fritz Knot Tier

L-19-534

\$420

Standard Version, 5mm, 29cm
Standard Version, 5mm, 40cm
Easy filling, suture size up to 1 metric
Easy handling, easy to tighten



New Dr Fritz Injection and Puncture Cannula

Interchangable needle (supplied 1mm x 26mm) with protection outer tube and inside tube with Luer Lock connector.

\$360

Small Animal: L-33-210 5mm, Working Length 33cm
Equine: L-47-210 5mm, Working Length 47cm
Spare Needle L-210 1mm x 26mm (2 Pack)

\$80



New Dr Fritz Equine 10mm Hand Instruments

Instruments are easily dismantable and rotatable through 360 Degrees. Each Instrument is supplied complete with jaw insert, sheath and ratchet handle.



131.11H05V

\$990

Mayo Scissors, Curved Left, Insulated
10mm, Working Length 37cm



131.48H21V

\$990

Dissecting Forceps Overholt, Turned, Insulated
10mm, Working Length 37cm



141.26H21V

\$990

Grasping Forceps Babcock,
10mm, Working Length 37cm



131.25H21V

\$990

Grasping Forceps Universal,
Double Action Jaws, Insulated
10mm, Working Length 37cm



131.49H21V

\$990

Grasping Forceps Atraumatic, Insulated
10mm, Working Length 37cm



131.32H21V

\$990

Grasping Forceps Claw, Double Action Jaws, Metal
10mm, Working Length 37cm






New Olympus WA29042B AI Endoscope with Sterilising Tray

Laparoscope/Endoscope supplied with Stopcock with Two Irrigation Ports, Instrument Tray and Sealing Caps. Angled ocular eyepiece. Ideal for Reproduction, ENT and Urology. Fully autoclavable.

\$9,460

Specifications:

- Working Length: 430mm
- Diameter: 8.6/9.8 Fr
- Instrument Channel: 6.4 Fr
- Direction of View: 7 Degree



New Standard Line 11mm Operating Laparoscope

Ideal for single Port Laparoscopic Procedures. Cannula and Trocar optional extra.

\$4,900

Specifications:

- Working Length: 312mm
- Diameter: 11mm
- Instrument Channel: 6.1mm



Olympus/Stryker 5mm Laparoscope

Ideal for Endoscopic Artificial Insemination in sheep and Small Animal/Equine Laparoscopic applications. Cannula and Trocar optional extra.

\$2,000

Specifications:

- Working Length: 300mm
- Diameter: 5mm
- Direction of View: 0 Degree or 30 Degree Available



Olympus/Stryker 10mm Laparoscope

Ideal for Endoscopic Artificial Insemination in sheep and Equine Laparoscopic applications and equine dentistry. Cannula and Trocar optional extra.

\$1,600

Specifications:

- Working Length: 300mm
- Diameter: 10mm
- Direction of View: 0 Degree or 30 Degree Available



Olympus/Stryker 10mm Extra Long Laparoscope

Ideal for Equine Laparoscopy or Dentistry. Cannula and Trocar optional extra.

\$3,000

Specifications:

- Working Length: 450mm
- Diameter: 10mm
- Direction of View: 0 Degree or 30 Degree Available

AESCULAP® CAIMAN®

ADVANCED BIPOLAR SEAL AND CUT TECHNOLOGY


A multipurpose vessel sealing instrument for Laparoscopic and Open surgery.

- Unique, tip first closure with patented floating hinge for vessel encapsulation
- Seals vessels up to 7mm Ø
- 26.5mm seal length, 23.5mm cut length
- 5mm Ø (Lap) and 12mm Ø (Open) in 24cm, 36cm and 44cm (articulating)

New, compact Seal and Cut generator with warranty and a readily available supply of hand pieces



Caiman GN200 Generator \$8,400



For more information download the B. Braun Caiman app through the Apple Store



| | | | | | |
|---------|--|------|--------|-----------------|---------|
| PL738SU | CAIMAN SEALING AND CUTTING INSTRUMENT | 5mm | 24cm | NON ARTICULATED | \$640 |
| PL740SU | CAIMAN SEALING AND CUTTING INSTRUMENT | 5mm | 36cm | NON ARTICULATED | \$640 |
| PL742SU | CAIMAN SEALING AND CUTTING INSTRUMENT | 5mm | 44cm | NON ARTICULATED | \$640 |
| PL739SU | CAIMAN SEALING AND CUTTING INSTRUMENT | 5mm | 24cm | ARTICULATED | \$730 |
| PL741SU | CAIMAN SEALING AND CUTTING INSTRUMENT | 5mm | 36cm | ARTICULATED | \$730 |
| PL743SU | CAIMAN SEALING AND CUTTING INSTRUMENT | 5mm | 44cm | ARTICULATED | \$730 |
| PL730SU | CAIMAN SEALING AND CUTTING INSTRUMENT | 12mm | 24cm | ARTICULATED | \$790 |
| PL731SU | CAIMAN SEALING AND CUTTING INSTRUMENT | 12mm | 44cm | ARTICULATED | \$790 |
| PL774SU | CAIMAN MARYLAND SEALING AND CUTTING INSTRUMENT | 5mm | 12.5cm | NON ARTICULATED | \$635 |
| PL775SU | CAIMAN MARYLAND SEALING AND CUTTING INSTRUMENT | 5mm | 17cm | NON ARTICULATED | \$635 |
| PL770SU | CAIMAN MARYLAND SEALING AND CUTTING INSTRUMENT | 5mm | 36cm | NON ARTICULATED | \$635 |
| PL772SU | CAIMAN MARYLAND SEALING AND CUTTING INSTRUMENT | 5mm | 44cm | NON ARTICULATED | \$635 |
| PL771SU | CAIMAN MARYLAND SEALING AND CUTTING INSTRUMENT | 5mm | 36cm | ARTICULATED | \$1,130 |
| PL773SU | CAIMAN MARYLAND SEALING AND CUTTING INSTRUMENT | 5mm | 44cm | ARTICULATED | \$1,130 |



New B. Braun Trocars & Cannulas

B | BRAUN

Ø 5mm

Sealing Units

Reusable

**EK080P \$380**Sealing unit
for trocars
Ø 5mm
(20 units)

Single unit price \$20

Trocar Pins**EK051R \$175**

triangular

EK045R \$175

conical sharp

EK052R \$175






triangular

EK046R \$175

conical sharp

Trocar Sleeves

Reusable

**EK010R \$265**smooth,
with
stopcock
60 mm**EK011R \$210**smooth,
without
stopcock
60 mm**EK012R \$320**threaded,
with
stopcock
60 mm**EK013R \$225**threaded,
without
stopcock
60 mm**EK014R \$265**smooth,
with
stopcock
110 mm**EK015R \$210**smooth,
without
stopcock
110 mm**EK016R \$320**threaded,
with
stopcock
110 mm**EK017R \$225**threaded,
without
stopcock
110 mm




Ø 10mm

Sealing Units

Reusable



EK083P \$140
Sealing unit
for trocars
Ø 10/12mm with
5 mm flap converter



EK086P \$90
Sealing unit
for trocars
Ø 10/12mm

Trocar Pins



EK062R \$175
triangular

EK056R \$175
conical sharp

EK063R \$175
triangular

EK057R \$175
conical sharp

Trocar Sleeves

Reusable



EK024R \$255
smooth,
with
stopcock
110 mm



EK025R \$210
smooth,
without
stopcock
110 mm



EK026R \$290
threaded,
with
stopcock
110 mm



EK027R \$225
threaded,
without
stopcock
110 mm

EK028R \$265
smooth,
with
stopcock
150 mm

B | BRAUN



Ø 12mm

Sealing Units

Reusable

**EK083P \$140**

Sealing unit
for trocars
Ø 10/12mm with
5 mm flap converter

**EK086P \$90**

Sealing unit
for trocars
Ø 10/12mm

Trocars Pins



110 mm

EK072R \$190
triangular

EK066R \$190
conical sharp

150 mm

EK073R \$190
triangular

EK067R \$190
conical sharp

Trocars Sleeves

Reusable



110 mm

EK034R \$265

smooth,
with
stopcock
110 mm

150 mm

EK038R \$300

smooth,
with
stopcock
150 mm



110 mm

EK035R \$210


smooth,
without
stopcock
110 mm



110 mm

EK036R \$310

threaded,
with
stopcock
110 mm



110 mm

EK037R \$220

threaded,
without
stopcock
110 mm

Sealing Units / Consumables

Reusable



EKO83P \$140
Sealing unit
for trocars
Ø 10/12mm with
5 mm flap converter



EKO87P \$410
Reducing converter
for trocars
Ø 10/12mm with
5 mm flap converter
(5 units)
Single unit price \$80



EKO86P \$90
Sealing unit
for trocars
Ø 10/12mm



EKO85P \$565
Sealing cap
for trocars
Ø 10/12mm
(20 units)
Single unit price \$30



EKO84P \$210
Cross slit valve
for trocars
Ø 10/12mm
(20 units)
Single unit price \$10

Trocar for mini-laparotomy (Hasson technique)

Reusable



EKO24R \$245
smooth,
with stopcock
110 mm



EKO64R \$175
blunt
(Hasson)
110 mm



EKO98R \$375
conical
attachment
for 10 mm
trocar sleeve

ARTG
Listed



Extra Long Cabinet

\$13,950



Andersen Anprolene Ethylene Oxide (EO) Sterilisation Cabinet

The most gentle, efficient and effective low temperature steriliser available. The Anprolene system offers unique plastic sterilisation bags and individual gas ampoules.

Features:

- **SAFE MICRO-DOSE** gas ampoule system – each sterilisation load uses less than 20gr of EO
- Unmatched sterilisation flexibility to sterilise a wide range of delicate **ENDOSCOPY** and **VETERINARY** instruments

Ideal For:

- Administration Sets
- Arthroscopes & Laparoscopes
- Biopsy & Grasping Forceps
- Bottles
- Catheters
- Cautery Cords & Tips
- Clothes
- Defibrillator Leads & Paddles
- Dermal Punches
- Drill Pieces
- Endotracheal Tubes
- Flexible Endoscopes
- Instruments
- Light Guide Cables
- Gloves
- Needles
- Ophthalmic Instruments
- Oscopes
- Respirators
- Surgical Staplers
- Syringes
- Tubing
- Video Endoscopy Cameras

Andersen Large Ampoule and Bag Refill Kit



In the Andersen system, the EO is supplied in a glass ampoule which is held within a break shield. Gas is delivered to the heart of the load, achieving the highest level of sterilisation for any given amount of gas.

The gas is released by breaking the neck of the ampoule. The EO boils, giving a very distinctive bubbling and cooling feeling as it does so. This provides very clear, positive feedback that the ampoule has been properly activated.



The packaged items, ready to be sterilised, are placed inside the supplied sterilisation bag. The cabinet holds the space around the sterilisation bag below room pressure. This ensures that any leak is from the room into the cabinet. Any EO released through the bag is taken to the outside through the exhaust tube at the back of the cabinet.

The bag volumes vary in size, giving great flexibility in the volume, size and shape of items that can be sterilised.

Large Ampoule and Bag Refill Kit

AN7916

\$850

consists of:

- 14 x 18cc Ampoules
- 14 x 40 Litre Sterilisation Bags
- Probe Strap

Andersen Dosimeter Chemical Integrator



Time, temperature and EO concentration are essential to successful sterilisation. Dosimeter integrators present visual assurance that all these parameters have been fulfilled. Dosimeters change colour in proportion to the dose of sterilant they contact. View results immediately after sterilisation. No laboratory culture is necessary.

Use one Dosimeter per cycle.

Features:

- Read results immediately
- Reliable
- Easy to use for every sterilisation load

AN87

\$115

Box of 25

Contact Us for Reference Sites



Austvet Endoscopy is committed to:

- Providing the highest level of Endoscopic service and repair
- Providing the fastest turnaround time
- Maximising our customers' satisfaction and productivity
- Providing training and technical support
- **NEW General Instrument Repair Service**
including tightening, realigning, wire brushing, sharpening, polishing, repinning, insulation testing, soldering and modifications to general and endoscopic instruments

OLYMPUS®

AOHUA



Don't Take Risks With Your Precious Endoscopy Equipment!



Endoscope Cleaning Brushes

Various sized brushes

\$10

- 1.0-1.5mm Instrument Channel/240cm Working Length
- 2.0-2.8mm Instrument Channel/220cm Working Length
- 2.0-2.8mm Instrument Channel/350cm Working Length



Matrix Wipes

Unique biofilm removing detergent wipes for pre-cleaning endoscopes and other medical devices. Ready to use and conveniently packaged wipes. Contain MatrixT biofilm removing solution (a patented formulation of surfactants specifically designed for removing biological soils and biofilms). Apertured wipes to ensure sufficient density for removing organic soiling.

\$10



Matrix Biofilm Remover Cleaning Solution

Matrix is a breakthrough in validated cleaning protocols. More effective than enzymatic cleaners, Matrix is the only proven solution designed to remove bacterial biofilms in which both bacteria and other microorganisms may be imbedded on flexible and rigid endoscopes and accessories which have accumulated over time due to insufficient cleaning.

\$200

5 Litre

Impacts 100% on Viable Bacteria to Provide Absolute Assurance That All Your Endoscopes and Instruments are Optimally Cleaned.



Medizyme Multi-Enzyme Cleaning Solution

Medizyme **Multi-Enzyme Cleaner** is totally compatible with all flexible and rigid endoscopes due to its neutral PH. Provides safe and effective cleaning against proteins, carbohydrates, lipids and mucopolysaccharides.

\$20

500ml

Also provides a combination of better detergency and good rinsing ability to ensure the best possible results.

\$110

5 Litre



OPAL Disinfectant

OPAL™ is an Ortho-Phthalaldehyde (OPA) based **high level liquid disinfectant** that features a reuse period of up to 21 days and a fast turnaround of 6-minute contact time at 20°C. Ready-to-use, no mixing or activation required. Effective against a wide array of micro-organisms. Suitable for use on all flexible and rigid endoscopes and accessories.

\$150

5 Litre

OPAL Test Strips

Opal® Test Strips are chemical indicators for use in determining whether the concentration of OPA, the active ingredient in OPAL®, is above or below the Minimum Recommended Concentration of 0.35%. It is recommended the OPAL® solution be tested before each disinfection cycle.

\$60

50 Test Strips

- Easy to use - 3 second dip time
- Results in 90 seconds
- Colour reference chart to confirm concentration amount. Large, easy to read indicator pad



Mr Clear Anti-Fog Solution

Convenient 6ml Anti-Fog solution which eliminates condensation at the distal end of rigid endoscopes/telescopes. Non-abrasive sponge allows for effective cleaning without distal lens damage.

\$15



Google or click on [Home | Olympus Academy](#) as a further resource for cleaning and maintenance of your endoscopes and all things endoscopic and consider booking your nursing staff into an **Olympus Reprocessing Workshop** in a capital city near you.



AUSTVET ENDOSCOPY



EO Sterilisation for the 21st Century

Austvet Endoscopy in conjunction with Steriflex are dedicated to the safe and effective sterilisation of delicate medical/veterinary devices. We are the sole Australian and New Zealand Distributors for the Andersen EO GAS Flexible Chamber Technology.

EO Gas not only offers the highest level of sterility but does so with negligible wear and tear to the instruments.

This not only saves money and time but provides you with peace of mind knowing that your equipment is sterile and that you are also helping the environment by reusing instruments and saving on landfill and waste.

EO (Ethylene Oxide) Gas is the only sterilisation modality that can cater for the sterilisation of the widest range of Endoscopes available today. **Refer to pages 32 & 33.**

We are also proud to represent

