

The World's Friendliest Pump

10k[®] Pump



Actual Size

■ Easy to Load ■ Quick to Prime ■ High Performance



 CONMED
LINVATEC

10k® Pump

The Simplest High-Performance Pump in Arthroscopy.

Ease of use and superior performance are the key characteristics that make the 10k Fluid System unique.

EASE OF USE

- Intuitive one-step/one-hand cassette loading
- Quick setup and priming
- Integrated pole mount works with any standard I.V. pole
- Compatible with video cart shelf or any flat surface

SUPERIOR PERFORMANCE

- Immediate response
- High flow
- Maintains distention under even the most demanding situations
- Robust autoclavable remote control

SPECIFICATIONS

- Maximum set pressure 150 mmHg
- Maximum flow rating 2000 ml/min
- Weight: 7 pounds
- Dimensions: 6.5" (16.5 cm) tall, 8" (20.3 cm) deep and 7.75" (19.7 cm) wide
- Voltage: 90 – 240 v

ORDER INFORMATION

Cat. No.	Description
10k	Console
10k100	Arthroscopy Inflow Tube Set
10k100D	Arthroscopy Inflow Tube with Colder Set
10k150	Arthroscopy Inflow/Outflow Tube Set
10k505	Laparoscopic Irrigation/Suction Tube Set
10k009	(Optional) Console Mounting Kit
C7115	Autoclavable Remote Control



10k Console



C7115 Autoclavable Remote Control




10k100 Arthroscopy Inflow Tube Set



10k100D Arthroscopy Inflow Tube with Colder Set



10k150 Arthroscopy Inflow/Outflow Tube Set



10k505 Laparoscopic Irrigation/Suction Tube Set



CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



11311 Concept Boulevard
Largo, FL 33773-4908
(727) 392-6464
Customer Service: 1-800-237-0169
FAX: (727) 399-5256
International FAX: +1 (727) 397-4540
email: customer_service@linvatec.com
www.linvatec.com