



AUSTVET ENDOSCOPY

PRODUCT CATALOGUE

2024

Est. 2004



## Endeavour for Better

**In-House Service & Repair.**

**Don't Take Risks With Your Precious Endoscopy Equipment!**

Please don't hesitate to call any time with any questions or queries  
or come and visit our showroom soon:

---

Unit 5 / 14 Lionel Rd, Mt Waverley VIC 3149

Maria: 0409 383 327 | Ennio: 0408 053 100 | Wally: 0438 084 418

[sales@austvetendoscopy.com.au](mailto:sales@austvetendoscopy.com.au) [www.austvetendoscopy.com.au](http://www.austvetendoscopy.com.au) [www.steriflex.com.au](http://www.steriflex.com.au)



Olympus Partnership	3
Olympus HD Video Endoscopy System	4
Endoscopic Auxiliaries	5
Olympus Flexible Video Endoscopes	6
Portable HD Video Endoscopy System	8
Portable Endoscope	9
Retrieval Devices	10
Therapeutic Devices	11
Fluid Administration & Therapy	12
Endoscopic Accessories	13
Arthroscopic Fluid Management	15
Otoscopy/Rhinology	16
Portable HD Video System	17
Arthroscopy	18
Laparoscopy	24
Caiman Vessel Sealing Device	27
Trocars & Cannulas	28
EO Gas Sterilisation	32
Cleaning Solutions	34
Service & Repair	35

**DISCLAIMER:** Please note that prices and availability are correct at the time of publishing but may be subject to change without notice.

**\*All prices exclude GST and freight**



**OLYMPUS®** True to Life

Olympus Australia's corporate philosophy is to make people and their much loved animal's lives healthier, safer and more fulfilled. One of their purposes is to *make the invisible, visible to see and treat disease*. They are on the cutting-edge of medical procedures and scientific discoveries and deliver a range of quality services and State of the Art endoscopy products.

The philosophy, "True to Life" sits in the heart of their beliefs and every day they develop and deliver medical technological solutions that elevate the standard of care to see and treat diseases with a full conviction that every life counts.

The medical instrument division encompasses a broad range of endoscopy and surgical instruments, with particular emphasis in the fields of:

Gastroscopy, colonoscopy, bronchoscopy and the entire spectrum of endoscopic and keyhole surgery, including laparoscopy, cystoscopy and arthroscopy.



**EVIS X1 Video System Centre CV-1500**  
Elevating the Standard of Endoscopy


**Austvet Endoscopy is the exclusive veterinary distributor of Olympus Endoscopy products in Australia and New Zealand.**

**Recognised as the world's leading supplier of medical endoscopy equipment, we are proud to offer the Olympus range of products to the veterinary industry and support their values and philosophies.**



Google or click on [continuum.olympusprofed.com](http://continuum.olympusprofed.com) as a further resource for cleaning and maintenance of your endoscopes and all things endoscopic and consider booking your nursing staff into an **Olympus Reprocessing Workshop** in a capital city near you.





## OLYMPUS®

The Olympus 190 Video System is cutting-edge technology for unsurpassable visualisation.

### Endoscopy Video System

#### Includes:

- Olympus CV-190 Video Processor with Keyboard and Scope Cable (MAJ-1430)
- Olympus CLV-190 Xenon Light Source with Water Bottle
- Endoscopy Trolley with Olympus OEV-261H High Definition 26" Monitor, Fully Adjustable Monitor Arm and Endoscope Hanger


#### Features:

- Full HD imaging capabilities when used in conjunction with H-Series video endoscopes.
- Narrow Band Imaging to enhance the visibility and characterisation of healthy and abnormal tissue and other structures on the mucosal surface
- Compatible with all Olympus Endoscopes & Surgical Camera Heads
- Image capture system is optional (see page 5)
- Video Camera head for Rigid Endoscopy optional

**\$29,950**



Camera Head Optional Extra





**New Endoscopy Trolley**

**\$1,950**


Trolley designed especially for Video Endoscopy.

**Features:**

- Heavy duty steel construction
- 4 tiers
- Medical grade castors with 100mm wheels
- Front castors with brakes
- Conceals all electrical wiring
- Power board 6-way, fully adjustable monitor arm and endoscope hanger

Dimensions: 550mm W x 550mm D x 1350mm H  
(Height excludes monitor)

**Monitor Optional Extra**



**New Microsoft Surface Laptop 4**

**\$3,490**

Microsoft Surface Laptop 4 with custom **HD** Image Capture Software.

**Features:**

- 13.5" display
- 2256 x 1504 screen resolution
- **HD** image and **HD** video capture capabilities
- Ability to input patient details (up to 80 characters)
- Supplied with SDI or DVI to USB converter

*Also note: Due to computer technology changing rapidly specifications and price may differ. Please enquire.*



## Specifications &amp; Compatibility Chart

Endoscope Model	Working Length	Diameter	Biopsy Channel	Angulation
URF-V2	670mm	2.8mm	1.2mm	2 Way
CYF-V2	380mm	5.4mm	2.2mm	2 Way
CYF-VH	380mm	5.5mm	2.2mm	2 Way
BF-P180	600mm	4.9mm	2mm	2 Way
BF-Q180	600mm	5.2mm	2mm	2 Way
BF-P190	600mm	4.2mm	2mm	2 Way
BF-Q190	600mm	4.8mm	2mm	2 Way
BF-H190	600mm	5.5mm	2mm	2 Way
GIF-XP150N	1100mm	5.5mm	2mm	4 Way
GIF-XP170N	1100mm	5.9mm	2mm	4 Way
GIF-XP180N	1100mm	5.5mm	2mm	4 Way
GIF-Q180	1030mm	8.8mm	2.8mm	4 Way
GIF-H180	1030mm	9.8mm	2.8mm	4 Way
GIF-H190	1030mm	9.2mm	2.8mm	4 Way
GIF-HQ190	1030mm	9.9mm	2.8mm	4 Way
PCF-H180AL	1680mm	11.8mm	3.2mm	4 Way
PCF-H190AL	1680mm	11.7mm	3.2mm	4 Way

*\*All Second Hand Endoscopes Pending Availability*



Image	Air / Water Facility	Processor CV-190	Processor CV-180	Price Second Hand
	No	✓	✗	\$11,900
	No	✓	✓	\$3,500
	No	✓	✓	\$6,900
	No	✓	✓	\$7,900
	No	✓	✓	\$7,900
	No	✓	✗	\$10,900
	No	✓	✗	\$10,900
	No	✓	✗	\$11,500
	Yes	✓	✓	\$9,500
	Yes	✓	✗	\$10,900
	Yes	✓	✓	\$13,900
	Yes	✓	✓	\$7,900
	Yes	✓	✓	\$10,900
	Yes	✓	✗	\$14,900
	Yes	✓	✗	\$14,900
	Yes	✓	✓	\$10,900
	Yes	✓	✗	\$14,900

OLYMPUS FLEXIBLE VIDEO ENDOSCOPES





Portable battery powered **HD** video-endoscopy system designed for ease of use for Veterinary surgeons performing upper airway and gastroscopic examinations.

Dimensions: 48cm x 42cm x 29cm

Weight: 17kg



### Aohua Vet OR1200 Li Endoscopy Video System

**\$16,500**

#### Includes:

- Battery Powered Video Processor with battery life indicator (over 3 hours of battery life)
- Microsoft Surface Book with 13.5" Screen
- Custom Made Carry Case
- Easy to use **HD** Image/Video Capture Software

### VET-9215 HD Video Endoscope

**\$16,500**

#### Specifications:

- Working Length: 1500mm
- Diameter: 9.2mm
- Biopsy Channel: 2.8mm
- Air/Water Facility
- 4 Way Angulation
- Includes Leak Tester, Water Bottle and Cleaning Accessories

### VET-9830 HD Video Endoscope

**\$16,900**

#### Specifications:

- Working Length: 3000mm
- Diameter: 9.8mm
- Biopsy Channel: 2.8mm
- Air/Water Facility
- 4 Way Angulation
- Includes Leak Tester, Water Bottle and Cleaning Accessories

### VET-6010 HD Video Endoscope

**\$16,900**

#### Specifications:

- Working Length: 1000mm
- Diameter: 6mm
- Biopsy Channel: 2mm
- Air/Water Facility
- 4 Way Angulation
- Includes Leak Tester, Water Bottle and Cleaning Accessories

### VET-3512 Video Endoscope

**\$15,900**

#### Specifications:

- Working Length: 600mm
- Diameter: 3.5mm
- Biopsy Channel: 1.2mm
- 2 Way Angulation
- Includes Leak Tester, Water Bottle and Cleaning Accessories





**AOHUA MVE-9215 Portable Endoscope**

**When portability is A MUST**

- Working Length: 150cm
- Diameter: 9.2mm
- Biopsy Channel: 2.8mm
- 4 Way Angulation
- LED Illumination

Scope with Small Monitor

**\$11,950**



**Clax Smart Fold Cart**

**\$490**

- Can be collapsed quickly and easily at the push of a button
- Has two large platforms, supplied standard with one collapsible storage box
- Has many useful features such as foot brake, holding device for boxes, upward folding platform and removable wheels for space-saving storage
- Made of injection moulded plastic and aluminium for long lasting use

**Specifications:**




Total capacity	60kg	Unit weight	7kg
Top load capacity	20kg	Bottom load capacity	40kg
Folding Handle type		Included	1 crate
Dimensions assembled	890 x 550 x 1030mm	Dimensions folded	670 x 470 x 110mm

**Great for transporting the Portable Equine Endoscopy System**



**New Oesophageal Rigid Foreign Body Removal Forceps**

Teeth face backwards to maximise grip and minimise damage to the oesophagus.

	Diameter at Jaws in mm	Length at Jaws in mm	Shaft Diameter in mm	Shaft Length in mm	Price
	13 x 9	34	7.5 x 6	640	<b>\$500</b>
	11.5 x 6	31	6.5 x 5	500	<b>\$500</b>
	5.6 x 4.2	15.5	4 x 3	400	<b>\$500</b>



**“Standard” Size Most Popular**


**New Flexible Retrieval Devices****Retrieval Grasping Basket**

For the extraction of round, slippery objects.

Paediatric: For 2.0mm > Instrument Channels, 180cm Length Basket Opening 32mm

Adult: For 2.8mm > Instrument Channels, 230cm Length Basket Opening 30mm

**Special order** for finer instrument channels < 1.2mm

**\$350****\$500****Retrieval Tripod Prongs**

For the collection of irregular shaped objects.

Paediatric: For 2.0mm > Instrument Channels, 180cm Length

Adult: For 2.8mm > Instrument Channels, 230cm Length

**\$150****SnareMaster Oval Loop**

Designed with optimal stiffness to precisely capture foreign bodies. Handle shaft is marked at 5mm increments to indicate the snare throw opening.

Paediatric: For 4.2mm > instrument Channels, 180cm Length Snare Opening 25mm

Adult: For 2.8mm > Instrument Channels, 230cm Length Snare Opening 30mm

For 2.8mm > Instrument Channels, 230cm Length Snare Opening 125mm (Jumbo)

For 2.8mm > Instrument Channels, 350cm Length Snare Opening 25mm


**\$50****A Cord for Hot Snare**

The cable connection between the Snare and Electrocautery which will enable cutting and coagulation of the tissue.

**\$100****Spider-Net Retrieval Device**

For the collection of round, blunt or otherwise hard to grasp foreign bodies. The “spider” net delivers a flexible, durable pouch eliminating the risk of accidental loss during retrieval.

Adult: For 2.8mm > Instrument Channels, 230cm Length Net Opening 3x6cm

**\$150****Ratigator Reusable Grasping Forceps**

Combines a rat tooth with alligator jaws into one instrument. 15mm jaw opening. Reusable & autoclavable.

Paediatric: For 2.0mm > Instrument Channels, 210cm Length


Adult: For 2.8mm > Instrument Channels, 230cm Length

For 2.8mm > Instrument Channels, 350cm Length

**\$700****Foreign Body Removal Hood**


This hood features a soft bell shaped hood for the removal of sharp or pointed gastric foreign bodies and folds backwards on insertion. The hood prevents exposure of the oesophageal and pharyngeal wall to injuries from the edges of the foreign body. Latex Free.

**\$100**



### New Flexible Biopsy Forceps


	Instrument Channel	Working Length	Price
<b>Serrated Cup Paediatric</b>	2.0mm	160cm	<b>\$50</b>
<b>Oval Cup Paediatric</b>	2.0mm	160cm	<b>\$20</b>
<b>Serrated Cup w Spike Paediatric</b>	2.0mm	160cm	<b>\$50</b>
<b>Serrated Cup Adult</b>	2.8mm	230cm	<b>\$20</b>
<b>Serrated Cup w Spike Adult</b>	2.8mm	230cm	<b>\$20</b>
<b>Oval Cup w Spike Adult</b>	2.8mm	230cm	<b>\$20</b>
<b>Oval Cup Adult</b>	2.8mm	230cm	<b>\$20</b>
<b>Oval Cup Adult</b>	2.8mm	400cm	<b>\$100</b>



### Flexible Cytology Brushes

Designed to enhance the efficiency of bronchoscopic and gastrointestinal examinations and improve collection of cytology samples. Features optimal brush hardness for better sample collection and thinner bristles for improved cell collection.

	Instrument Channel	Working Length	Brush Size	Price
<b>Bronchial</b>	1.2mm	115cm	2mm	<b>\$50</b>
	2.0mm	120cm	3mm	<b>\$50</b>
	2.0mm	120cm	2mm	<b>\$50</b>
<b>Gastrointestinal</b>	2.0mm	200cm	3mm	<b>\$50</b>
	2.8mm	200cm	3mm	<b>\$50</b>



### Injector Needle

Designed to provide proper flexibility and rigidity to ensure safe piercing and excellent puncture performance. Ideal for injecting Formalin to treat Ethmoid Haematomas in horses. Assorted needle lengths and gauges. Large range of needle gauges and lengths for paediatric and adult scopes, please enquire

Instrument Channel	Working Length	Needle Diameter	Needle Length	Price
2.0mm	200cm	23G	5mm	<b>\$80</b>
2.8mm	230cm	23G	6mm	<b>\$80</b>
2.8mm	240cm	21G	8mm	<b>\$80</b>



### Olympus EZDilate Inflation Device

**\$150**

Used while performing balloon dilation procedures to inflate the balloon, monitor the pressure within the balloon and deflate the balloon. Easy to read luminescent dial and quick release lock lever mechanism. Maximum pressure of 30 atmospheres.




### Olympus EZDilate Multi-Stage Endoscopic Balloon Dilators

**\$350**

Designed for oesophageal dilation under direct endoscopic vision. Low compliance balloon material produces dependable, controllable radial dilation. Low profile balloon design assures smooth accessibility into difficult strictures.

Assorted Sizes: For 2.8mm > Instrument Channels, from 6mm to 20mm Diameter

Inflation diameter (mm)	6, 7, 8	8.5, 9.5, 10.5	11, 12, 13	13.5, 14.5, 15.5	16, 17, 18	18, 19, 20
<b>Inflation pressure (atm)</b>	2, 5.5, 9.5	2, 5.5, 9.5	2, 3.5, 5.5	2.5, 4, 6	2.5, 4, 5.5	2.5, 3.5, 5
<b>Balloon length (mm)</b>	65	65	55	55	55	55



### Mila Double/Triple Lumen Flushing Catheter

Flushing catheter for guttural pouch lavage. Can also be used as a sampling catheter for sterile specimen retrieval via an endoscope. The set contains a glycol plug in the outer catheter and an inner catheter for retrieval of the sterile specimen.

**1.8mm Diameter/250cm Length/Double Lumen**

**\$100**

**2.3mm Diameter/190cm Length/Triple Lumen**

**\$100**



### Equine Uterine Flushing Catheter

Used for flushing and irrigation, administration of medications, embryo collection and transfer. Autoclavable. Sterile.

**12mm Diameter/80cm or 150cm Length**

**\$150**



### Endoscopic BAL Small Animal/Equine Single Lumen Flushing Catheter

An endoscopically guided catheter for lung lavage or introduction into the guttural pouch to flush out pouch debris. Visual indicators at 2.5cm, 5cm and 7.5cm from the distal tip to aid in introduction. Sterile.

**1.8mm Diameter/120cm Length**

**\$50**

**1.8mm Diameter/250cm Length**

**\$50**

**2.3mm Diameter/220cm Length**

**\$15**

**2.3mm Diameter/400cm Length**

**\$20**



### Silicone Stomach Tube

Three strategically placed fenestrations at the distal tip provide improved safety and effectiveness. Includes syringe and funnel adaptor. Autoclavable. Non-Sterile. (Internal dimensions listed).

**4.5mm/91cm Length or 5.0mm/91cm Length**

**\$80**

**5.0mm/183cm Length or 6.5mm/150cm Length**

**\$80**

**6.5mm/300cm Length or 9.5mm/300cm Length**

**\$120**

**11mm/300cm Length or 16mm/305cm Length**

**\$150**

**10mm/300cm Length (Standard Line)**

**\$100**



### Broncho-Alveolar Lavage Catheter (Silicone) BAL TUBE

The BAL catheter is used for lung lavage of selected segments of the lung fields in horses. Equipped with a distal balloon that can be inflated to prevent retrograde fluid flow and a proximal two-way fitting allows simultaneous connection for both irrigation and aspiration. Sterile and Autoclavable

(The plastic proximal two-way fitting and three-way stopcock are not autoclavable).

**10mm Diameter/240cm Length**

**\$180**

**10mm Diameter/300cm Length**

**\$200**

**Three-Way Replacement Stopcock**

**\$20**





**Standard Single Scope Hanger**

Requires only two screws for mounting.  
Dimensions: 145mm W x 135mm D

**\$150**



**Olympus, Aohua and Pentax Water Bottles**

Pentax Water Bottle to suit all Pentax Endoscopes  
Olympus Water Bottle Autoclavable for 140/160/180 Series  
Olympus Endoscopes  
Olympus Water Bottle Autoclavable for 100 Series  
Olympus Endoscopes

**\$500**

**\$300**

**\$200**

**Second Hand Bottles Also Sold Pending Availability**



**Olympus, Aohua and/or Pentax Leakage Testers**

Designed to detect and identify leaks from your flexible endoscopes. Stops potential fluid invasion thereby preventing the possibility of costly repairs, cross infection and downtime.

**\$300**

**ALL FLEXIBLE ENDOSCOPES MUST BE LEAK TESTED PRIOR TO EVERY SUBMERSION TO ENSURE YOU, YOUR PATIENTS AND YOUR SCOPES ARE PROTECTED.**



**Askir Endoscopic Suction Pump & Accessories**

Compact, high performance unit designed for both endoscopic aspiration and surgical suction. Unit is supplied complete with polycarbonate autoclavable 1 litre jar with overflow valve and antibacterial filter. Pump is equipped with aspiration regulator and vacuum indicator. Designed for almost continuous use up to 120 minutes. Maximum suction flow 30 litres per minute. Portable and lightweight.

**\$690**

**Accessories:**



Mucus Specimen Trap  
option of 40cc and 70cc

**\$15**



Suction Handpiece Frazier  
2.6mm, 3.3mm or 4.0mm

**\$30**



Suction Handpiece Yankauer 29cm  
Standard or Paediatric

**\$40**



### NEW AE-L100 Endoscopic LED Light Source (No Air Pump)

**\$1,390**

Ideal for use with all rigid endoscopes including Arthroscopes, Laparoscopes, Oscopes, Rhinoscopes and Cystoscopes.

#### Features:

- Extremely compact (30cm x 20cm x 10cm)
- Expected lifetime of LED is over 40,000 hours
- Extremely bright light
- Supplied with light guide cable



### For Rigid Endoscopes

### LED Portable High-Intensity Light Handles for Rigid and Flexible Endoscopes Now with Variable Brightness

LED Portable Light Handles are an easy option for quick examinations in the consult room or for fieldwork. They feature high-intensity illumination and are a great light source backup. Handles attach to most rigid and flexible endoscopes.



### For Flexible Endoscopes

#### Features:

- Lithium-Ion rechargeable battery
- Extremely compact
- High-intensity white light with variable brightness

#### Kit Includes:

- LED light source handle
- Rechargeable Lithium-Ion battery
- Battery charger

LED Light Handle for Rigid Endoscopes **\$250**

LED Light Handle for Flexible Endoscopes **\$300**



### New Sterile Camera Sleeve with Telescope Adaptor

Disposable surgical sterile video camera sleeve/cover with enclosed adaptor for telescope connection. Safe, quick and easy application. Latex free. Offers 100% sterility without chemicals to reduce risk of infection. Excellent equipment protection. Dimensions: 14cm W x 285cm L

**\$35** Each  
or  
**\$1,200**  
Box of 40



### Olympus UHI-3 35 Litre Insufflator

**\$795**

Olympus insufflator which features preset pressure and flow rate, leakage warning display and overpressure warning display and alarm. Supplied with hose for gas bottle. Designed to support a wide range of Laparoscopic insufflation applications.



### New Laparoscopic Insufflator Tubing

**\$15**

Disposable CO2 Insufflator Filter with 3m Tubing. Has universal fitting however any additional connectors required can be supplied at no charge. Filtration area 24cm<sup>2</sup> for added protection and greater gas flow. Sterile and Latex free.



### New Dr Fritz Veress Needles

Veress Needles equipped with a spring loaded obturator to insufflate the abdomen in  
 100mm **L-19-100** **\$220**  
 150mm **L-19-150** **\$300**  
 Laparoscopic surgery.



**New Linvatec 10k Pump. The Simplest High-Performance Pump in Arthroscopy**

**Ease of Use**

- Intuitive one-step/one-hand cassette loading
- Quick Setup and priming
- Integrated pole mount works with any standard I.V. pole
- Compatible with video cart shelf or any flat surface

**Superior performance**

- Immediate response
- High flow
- Maintains distension under even the most demanding situations
- Robust autoclavable remote control

**Specifications**

- Maximum set pressure 150mmHg
- Maximum flow rating 2000ml/min
- Weight: 3kg
- Dimensions: 16.5cm tall, 20.3cm deep and 19.7cm wide
- Voltage 90–240v



**10k100 Endoscopy Tube Set**

**\$700**  
box of 10


**\$5,500**



**New Dr Fritz Pressure Infusion Cuff**

500ml - 1 Litre (with manometer)	<b>AR-08-100N</b>	<b>\$270</b>
3 - 5 Litre (with manometer)	<b>AR-08-300N</b>	<b>\$770</b>





### New Standard Veterinary Rhinoscope Set

\$5,900

Rhinoscope set supplied with two sheaths, with the larger sheath allowing both biopsy and irrigation.

Features universal eyepiece for compatibility with all Endoscopic video cameras.

Multi-purpose rigid endoscope which can also be used for:

**Arthroscopy**

**Avian Endoscopy**

**Cystoscopy**

**Laparoscopy**

**Thoracoscopy**

**Vaginoscopy**

#### System Includes:

- Rhinoscope 2.7mm 30 Deg (190mm Working Length)
- 3.7mm/11 Fr Sheath (no biopsy facility – only irrigation)
- 4.7mm/14 Fr Sheath (biopsy & irrigation)
- High-intensity LED Portable Light Handle (with variable brightness)
- Cleaning Brush & Custom Carry Case
- Supplied with your choice of 1.6mm Biopsy or Grasping Forceps (see page 17)

Otoscope not included (please refer to below)

### New Standard Veterinary Otoscope Set

\$2,100

Complete, all-in-one veterinary otoscope. Ideal for quick examinations in the consult room.

Large field of view and depth of field.

Also used for foreign body removal and/or biopsy and quick diagnostics.

#### Set Includes:

- Otoscope 5.5mm 0 Deg (100mm Working Length)
- High-intensity LED Portable Light Handle (with variable brightness)
- Two Way Stopcock (Removable)
- Cleaning Brush & Carry Case (refer to image of carry case on page 16)
- Supplied with your choice of 1.6mm Biopsy or Grasping Forceps (see page 17)



**Essential Instrument for Every Clinic**





**New All-in-One Portable Video System (AE-1080 HD)**

**\$7,950**

This all-in-one portable **HD** Video System is a complete, user friendly solution for all types of rigid endoscopy for the practice as well as the mobile vet.

**Features:**

- 15.5" full **HD** display
- 3 metre long camera and light lead cables
- Image/video capture facilitated from camera head
- Multiple video outputs for slave monitors

**System Includes:**

- 15.5" full **HD** display
- Video processor
- Camera head
- LED light source
- Light Lead
- Custom carry case

**New 1.6mm/5fr Semi Flexible Biopsy/Retrieval Forceps**

**\$390**


Will suit rigid ENT/Otoscopes with instrument channels 1.6mm and above.

Working length: 270mm.




**48404**

Ultra slim 1.6mm/5fr semi flexible biopsy forcep. Ideal for biopsy of ear, nose or throat.




**48405**

Ultra slim 1.6mm/5fr semi flexible grasping forcep. Ideal for removal of foreign bodies from ear canals and nasal cavities.



**48405-B**

Ultra slim 1.6mm/5fr semi flexible dual purpose biopsy/retrieval forcep.



**New Olympus A37026A ENT/Cystoscopy Endoscope & Accessories**

**\$9,520**

Compact slim design and lightweight, this scope is the ideal choice for ENT and Cystoscopy. Integrated stopcock saves time and simplifies procedure. Atraumatic smooth beak ensures optimum patient comfort.

Supplied with instrument tray and sealing caps.

**Specifications:**

- 7 Degree field of View
- 2.6mm Diameter
- 160mm Working Length
- 1.4mm Instrument Channel
- 45 Degree Angled Ocular
- Fully Autoclavable

**Accessories:**

**New Flexible Grasping Forceps**  
1mm / 280mm Working Length

**\$700**

**New Flexible Biopsy Forceps**  
1mm / 280mm Working Length


**\$600**






2.7mm

**A70963A** Telescope,  
**\$3,300** Autoclavable,  
 30° direction of view,  
 Wide angle field of view,  
 70mm working length



**WA70992A** Trocar tube,  
**\$1,200** Rotatable,  
 1 Stopcock,  
 Standard flow,  
 for 70mm telescopes




**A70978A** Trocar spike,  
**\$290** triangular,  
 for 70mm telescopes



**A70979A** Trocar spike,  
**\$290** blunt,  
 for 70mm telescopes




**A00507A** Sealing, set  
**\$80** Set of spare sealing rings:  
 Contains 2 large sealings  
 (Transparent) and 5 small sealings (black)




**A0460** Light-guide adapter,  
**\$80** for OES 4000 light-guide cables

**00332** Light-guide adapter,  
**\$30** for Wolf light-guide cables


**A5333** Eyepiece cup,  
**\$110** for telescopes




4mm



**A70941A** Telescope,  
**\$3,600** Autoclavable,  
 30° direction of view,  
 115° field of view,  
 160mm working length



**WA70952A** Trocar tube,  
**\$1,200** Rotatable,  
 1 Stopcock,  
 Standard flow




**A70954A** Trocar spike,  
**\$290** Triangular



**WA70953A** Trocar tube,  
**\$1,200** Rotatable,  
 2 Stopcocks,  
 Standard flow



**A70955A** Trocar spike,  
**\$290** Blunt




**New Standard Line 3mm ENT/Arthroscope**

Arthroscope supplied with 4mm Cannula and Sharp and Blunt Trocar.

**\$1,750**

**Specifications:**

- Working Length: 175mm (with cannula 135mm)
- Diameter: 3mm
- Direction of View: 30 Degree



**New Standard Line 4mm ENT/Arthroscope**

Arthroscope supplied with 5mm Cannula and Sharp and Blunt Trocar.

**\$1,750**

**Specifications:**

- Working Length: 175mm (with cannula 135mm)
- Diameter: 4mm
- Direction of View: 30 Degree





### New Dr Fritz Small Animal Arthroscopy Hand Instruments

Instruments provide easy handling and an ergonomic design with handles. All instruments are specifically designed for Small Animal Arthroscopy.



#### AR-17-03S

Bayonet Knife with Handle, 2.7mm

\$400



#### AR-17-05S

Banana Knife with Handle, 2.7mm

\$400



#### AR-17-12S

Mikro Picking Knife with Handle, 2.7mm

\$400



#### AR-17-19S

Hook Knife with Handle, 2.7mm

\$400



#### AR-17-06S

Curette Curved with Handle, 2.7mm

\$400



#### AR-17-07S

Ring Curette with Handle, 2.7mm

\$400



#### AR-17-11S

Meniscutome Smilie with Handle, 2.7mm

\$400



#### AR-17-09S

Hook Probe with Handle, 2.2mm

\$400



#### AR-17-13

Special Elevator with Handle, 2.7mm (Metal)

\$400



#### I-17-822 (SET)

Milling Drill  
2.2mm, Working Length 12cm,  
Rotating Handpiece and Long Lasting  
Drill Steel (Dismantable)

\$750



Replacement Drill Steel (for Milling Drill)

\$300



#### S110-1935

Protection & Transportation Sheath for  
Small Arthroscopes 12cm, Autoclavable

\$50

**New Dr Fritz Small Animal Arthroscopy  
Punches, Forceps & Scissors**



**AS-17-521** **\$1,580**  
Mini Scissors for Arthroscopy, 2mm,  
Working Length 11cm



**AT-19-231** **\$2,150**  
Mini Rongeur & Biopsy Forceps, 2.2/2.7mm,  
Working Length 11cm



**AS-17-632** **\$1,580**  
Universal Rongeur & Grasping Forceps, 2.7mm,  
Working Length 12cm



**AS-17-662** **\$1,950**  
Alligator Grasping Forceps with Ratchet, 2.7mm,  
Working Length 12cm



**AS-17-621** **\$1,580**  
Arthroscopic Hook Scissor, 2.7mm,  
Working Length 12cm




**AS-17-672** **\$1,760**  
Arthroscopic Punch, 2.7mm,  
Working Length 12cm



**AS-17-641** **\$1,580**  
Biopsy & Spoon Forceps, 2.7mm,  
Working Length 12cm

**New Dr Fritz Small Animal Working Cannula Set**

This system allows the specific exchange of different working cannulas. Beginning with the smallest diameter, the vet can always change to a bigger diameter by using the changing rods supplied. Cannulas are specifically designed for Small Animal Arthroscopy.



**I-17-755** **\$1,170**  
**Working Cannula Set**

Consists of:

- 3 Cannulas for Elbow 2.3mm, 2.9mm and 3.5mm, Working Length 3cm – 5.5cm
- 2 Cannulas for Shoulder 2.9mm and 3.5mm, Working Length 7cm – 8cm
- 2 Sharp Trocars
- 2 Changing Rods, 2.3mm and 2.8mm



### New Dr Fritz Equine Arthroscopy Hand Instruments

Instruments have been specifically designed for Equine Arthroscopy.



#### I-18-100

Osteotome Curved with Handle, 4mm

\$240



#### I-18-110

Osteotome Straight with Handle, 4mm

\$240



#### I-18-520

Hook Probe with Handle, 3mm, Working Length 12cm

\$310



#### I-18-521

Special Knife with Handle, 3mm, Break Safe Special Cut

\$390



#### I-18-601

Curette Straight with Handle, Oval Cup 6 x 4mm

\$300



#### I-18-621

Curette, 35 Degree Oblique Left, Oval Cup, 6 x 4mm

\$320



#### I-18-622

Curette, 35 Degree Oblique Right, Oval Cup, 6 x 4mm

\$320



#### I-18-518

Egress Cannula and Stopcock with Sharp Trocar, 2.7mm, Working Length 7cm  
(Blunt Trocar additional \$150)

\$320



#### I-18-519

Egress Cannula with Sharp Trocar and Stopcock, 4mm, Working Length 7cm

\$320



#### I-18-635

Special Egress Curette Straight, Oval Cup, 5.5 x 3mm and Leur Lock (f) Suction Connector, Working Length 12cm

\$550



#### I-18-536

Special Elevator for Coffin Joint, Working Length 8cm

\$500



#### I-18-517

Elevator Curved with Handle, 6mm

\$350



#### I-18-534

Special Elevator for Sesamoide Fractures, Curved, Long 4mm

\$400



**New Dr Fritz Equine Arthroscopy Line Instruments**

**“Standard Line” Rongeur**

**I-18-6611**

**\$1,580**

Straight, Fenestrated, Oval Cup 4 x 8mm, Working Length 11cm

**I-18-662**

**\$1,580**

Angled Upwards, Fenestrated, Oval Cup 4 x 8mm, Working Length 11cm



**“Power Line” Rongeur**

**I-18-973**

**\$1,800**

Straight, Heavy Durable Handle, Oval Cup 5 x 10mm, Working Length 15cm

**I-18-974**

**\$1,800**

Straight, Heavy Durable Handle, Fenestrated, Small Tip, Oval Cup 5 x 10mm, Working Length 15cm



**New Ferris Smith Equine Rongeurs**

**\$1,100**

Premium quality German made Rongeurs. Features straight and angled upward tip with 180mm working length. Lifetime warranty.

**Sizes Available:**


- 3mm W x 10mm L Tip (Straight & Angled Upward)
- 4mm W x 10mm L Tip (Straight & Angled Upward)
- 5mm W x 10mm L Tip (Straight & Angled Upward)
- 6mm W x 10mm L Tip (Straight only)





### New Dr Fritz Small Animal 5mm Hand Instruments

Instruments are easily dismantable and rotatable through 360 Degrees. Each Instrument is supplied complete with jaw insert, sheath and ratchet handle.



#### 531.61H01V

Biopsy Forceps, Double Action Jaws, With Spikes 5mm, Working Length 33cm

\$825

#### 541.61H01V

Biopsy Forceps, Double Action Jaws, With Spikes 5mm, Working Length 43cm

\$825



#### 531.01H05V


Mini Metzenbaum Scissors, Curved Left, Insulated 5mm, Working Length 33cm

\$825

#### 541.01H05V

Mini Metzenbaum Scissors, Curved Left, Insulated 5mm, Working Length 43cm

\$825



#### 531.22H11V

Grasping Forceps Universal, Insulated 5mm, Working Length 33cm

\$825

#### 541.22H11V

Grasping Forceps Universal, Insulated 5mm, Working Length 43cm

\$825



#### 531.40H01V


Dissecting Forceps Maryland, Turned, Insulated 5mm, Working Length 33cm

\$825

#### 541.40H01V

Dissecting Forceps Maryland, Turned, Insulated 5mm, Working Length 43cm

\$825



#### 531.30H12V


Grasping Forceps Claw, Straight, Metal 5mm, Working Length 33cm

\$825

#### 541.30H12V

Grasping Forceps Claw, Straight, Metal 5mm, Working Length 43cm

\$825



#### 531.27H11V

Grasping Forceps Babcock 5mm, Working Length 33cm

\$825





**New Dr Fritz Needle Holder with Handle**

Needle self positioning, axial handle.  
Small Animal:



**L-19-562V** **\$1,560**

Left Turned 5mm,  
Working Length 22cm

**L-19-563V** **\$1,620**

Angled Upwards 5mm,  
Working Length 22cm

**New Dr Fritz Knot Tier**



**L-19-534** **\$420**

Standard Version, 5mm, 29cm

**L-19-544** **\$450**

Standard Version, 5mm, 40cm  
Easy filling, suture size up to 1 metric  
Easy handling, easy to tighten



**New Dr Fritz Injection and Puncture Cannula**

Interchangeable needle (supplied 1mm x 26mm) with protection outer tube and inside tube with Luer Lock connector.

Small Animal: **L-33-210** 5mm, Working Length 33cm **\$450**

Equine: **L-47-210** 5mm, Working Length 47cm **\$620**

Spare Needle **L-210** 1mm x 26mm (2 Pack) **\$150**






### New Olympus WA29042B AI Endoscope with Sterilising Tray

Laparoscope/Endoscope supplied with Stopcock with Two Irrigation Ports, Instrument Tray and Sealing Caps. Angled ocular eyepiece. Ideal for Reproduction, ENT and Urology. Fully autoclavable.

**\$9,900**

#### Specifications:

- Working Length: 430mm
- Diameter: 8.6/9.8 Fr
- Instrument Channel: 6.4 Fr
- Direction of View: 7 Degree



### New Standard Line 11mm Operating Laparoscope

Ideal for single Port Laparoscopic Procedures. Cannula and Trocar optional extra.

**\$4,900**

#### Specifications:

- Working Length: 312mm
- Diameter: 11mm
- Instrument Channel: 6.1mm



### Olympus/Stryker 5mm Laparoscope

Ideal for Endoscopic Artificial Insemination in sheep and Small Animal/Equine Laparoscopic applications. Cannula and Trocar optional extra.

**\$2,000**

#### Specifications:

- Working Length: 300mm
- Diameter: 5mm
- Direction of View: 0 Degree or 30 Degree Available



### Olympus/Stryker 10mm Laparoscope

Ideal for Endoscopic Artificial Insemination in sheep and Equine Laparoscopic applications and equine dentistry. Cannula and Trocar optional extra.

**\$1,600**

#### Specifications:

- Working Length: 300mm
- Diameter: 10mm
- Direction of View: 0 Degree or 30 Degree Available



### Olympus/Stryker 10mm Extra Long Laparoscope

Ideal for Equine Laparoscopy or Dentistry. Cannula and Trocar optional extra.

**\$5,000**

#### Specifications:

- Working Length: 450mm
- Diameter: 10mm
- Direction of View: 0 Degree or 30 Degree Available

# AESCULAP® CAIMAN®


ADVANCED BIPOLAR SEAL AND CUT TECHNOLOGY

A multipurpose vessel sealing instrument for Laparoscopic and Open surgery.

- Unique, tip first closure with patented floating hinge for vessel encapsulation
- Seals vessels up to 7mm Ø
- 26.5mm seal length, 23.5mm cut length
- 5mm Ø (Lap) and 12mm Ø (Open) in 24cm, 36cm and 44cm (articulating)



**New, compact Seal and Cut generator with warranty and a readily available supply of hand pieces**



For more information download the B. Braun Caiman app through the Apple Store

**Caiman GN200 Generator \$8,600**



PL738SU	CAIMAN SEALING AND CUTTING INSTRUMENT	5mm	24cm	NON ARTICULATED	<b>\$670</b>
PL740SU	CAIMAN SEALING AND CUTTING INSTRUMENT	5mm	36cm	NON ARTICULATED	<b>\$670</b>
PL742SU	CAIMAN SEALING AND CUTTING INSTRUMENT	5mm	44cm	NON ARTICULATED	<b>\$670</b>
PL739SU	CAIMAN SEALING AND CUTTING INSTRUMENT	5mm	24cm	ARTICULATED	<b>\$770</b>
PL741SU	CAIMAN SEALING AND CUTTING INSTRUMENT	5mm	36cm	ARTICULATED	<b>\$770</b>
PL743SU	CAIMAN SEALING AND CUTTING INSTRUMENT	5mm	44cm	ARTICULATED	<b>\$770</b>
PL730SU	CAIMAN SEALING AND CUTTING INSTRUMENT	12mm	24cm	ARTICULATED	<b>\$830</b>
PL731SU	CAIMAN SEALING AND CUTTING INSTRUMENT	12mm	44cm	ARTICULATED	<b>\$830</b>
PL774SU	CAIMAN MARYLAND SEALING AND CUTTING INSTRUMENT	5mm	12.5cm	NON ARTICULATED	<b>\$660</b>
PL775SU	CAIMAN MARYLAND SEALING AND CUTTING INSTRUMENT	5mm	17cm	NON ARTICULATED	<b>\$660</b>
PL770SU	CAIMAN MARYLAND SEALING AND CUTTING INSTRUMENT	5mm	36cm	NON ARTICULATED	<b>\$660</b>
PL772SU	CAIMAN MARYLAND SEALING AND CUTTING INSTRUMENT	5mm	44cm	NON ARTICULATED	<b>\$660</b>
PL771SU	CAIMAN MARYLAND SEALING AND CUTTING INSTRUMENT	5mm	36cm	ARTICULATED	<b>\$1,180</b>
PL773SU	CAIMAN MARYLAND SEALING AND CUTTING INSTRUMENT	5mm	44cm	ARTICULATED	<b>\$1,180</b>





Ø 5mm

**Sealing Units**

Reusable




**EK080P \$400**

Sealing unit for trocars Ø 5mm (20 units)

Single unit price \$25

**Trocar Pins**



**EK051R \$185**  
triangular


**EK045R \$185**  
conical sharp

**EK052R \$185**  
triangular

**EK046R \$185**  
conical sharp

**Trocar Sleeves**

Reusable




**EK010R \$280**  
smooth, with stopcock  
60 mm



**EK011R \$225**  
smooth, without stopcock  
60 mm




**EK012R \$340**  
threaded, with stopcock  
60 mm




**EK013R \$240**  
threaded, without stopcock  
60 mm




**EK014R \$280**  
smooth, with stopcock  
110 mm



**EK015R \$225**  
smooth, without stopcock  
110 mm



**EK016R \$340**  
threaded, with stopcock  
110 mm



**EK017R \$240**  
threaded, without stopcock  
110 mm



Ø 10mm

**Sealing Units**

Reusable




**EK083P \$150**  
Sealing unit  
for trocars  
Ø 10/12mm with  
5 mm flap converter



**EK086P \$95**  
Sealing unit  
for trocars  
Ø 10/12mm

**Trocar Pins**



110 mm

**EK062R \$185**  
triangular

**EK056R \$185**  
conical sharp

150 mm

**EK063R \$185**  
triangular

**EK057R \$185**  
conical sharp

**Trocar Sleeves**

Reusable



**EK024R \$270**  
smooth,  
with  
stopcock  
110 mm




**EK025R \$225**  
smooth,  
without  
stopcock  
110 mm



**EK026R \$305**  
threaded,  
with  
stopcock  
110 mm



**EK027R \$240**  
threaded,  
without  
stopcock  
110 mm



**EK028R \$280**  
smooth,  
with  
stopcock  
150 mm





Ø 12mm

**Sealing Units**

Reusable




**EK083P \$150**  
Sealing unit for trocars Ø 10/12mm with 5 mm flap converter



**EK086P \$95**  
Sealing unit for trocars Ø 10/12mm

**Trocar Pins**



**EK072R \$200**  
triangular

**EK066R \$200**  
conical sharp

**EK073R \$200**  
triangular

**EK067R \$200**  
conical sharp

**Trocar Sleeves**

Reusable



**EK034R \$280**  
smooth, with stopcock 110 mm



**EK035R \$225**  
smooth, without stopcock 110 mm



**EK036R \$330**  
threaded, with stopcock 110 mm



**EK037R \$235**  
threaded, without stopcock 110 mm



**EK038R \$315**  
smooth, with stopcock 150 mm



Sealing Units / Consumables

Reusable



**EKO83P \$150**  
Sealing unit for trocars Ø 10/12mm with 5 mm flap converter



**EKO87P \$435**  
Reducing converter for trocars Ø 10/12mm with 5 mm flap converter (5 units)  
*Single unit price \$90*



**EKO86P \$95**  
Sealing unit for trocars Ø 10/12mm



**EKO85P \$600**  
Sealing cap for trocars Ø 10/12mm (20 units)  
*Single unit price \$35*




**EKO84P \$225**  
Cross slit valve for trocars Ø 10/12mm (20 units)  
*Single unit price \$15*

Trocar for mini-laparotomy (Hasson technique)

Reusable



**EKO24R \$260**  
smooth, with stopcock 110 mm



**EKO64R \$185**  
blunt (Hasson) 110 mm



**EKO98R \$395**  
conical attachment for 10 mm trocar sleeve

ARTG  
Listed



Extra Long Cabinet

**\$14,500**



**ANPRO**  
ANDERSEN  
PRODUCTS

 **Steriflex**  
SAFE STERILISATION SOLUTION



### Andersen Anprolene Ethylene Oxide (EO) Sterilisation Cabinet AN74IX

The most gentle, efficient and effective low temperature steriliser available. The Anprolene system offers unique plastic sterilisation bags and individual gas ampoules.

#### Features:

- **SAFE MICRO-DOSE** gas ampoule system – each sterilisation load uses less than 20gr of EO
- Unmatched sterilisation flexibility to sterilise a wide range of delicate **ENDOSCOPY and VETERINARY** instruments

#### Ideal For:

- Administration Sets
- Arthroscopes & Laparoscopes
- Biopsy & Grasping Forceps
- Bottles
- Catheters
- Cautery Cords & Tips
- Clothes
- Defibrillator Leads & Paddles
- Dermal Punches
- Drill Pieces
- Endotracheal Tubes
- Flexible Endoscopes
- Instruments
- Light Guide Cables
- Gloves
- Needles
- Ophthalmic Instruments
- Oscopes
- Respirators
- Surgical Staplers
- Syringes
- Tubing
- Video Endoscopy Cameras
- Anything and Everything Except Food & Drugs

**PROTECT YOUR PATIENTS, INSTRUMENTS AND ENVIRONMENT**



**Andersen Large Ampoule and Bag Refill Kit**



In the Andersen system, the EO is supplied in a glass ampoule which is held within a break shield. Gas is delivered to the heart of the load, achieving the highest level of sterilisation for any given amount of gas.

The gas is released by breaking the neck of the ampoule. The EO boils, giving a very distinctive bubbling and cooling feeling as it does so. This provides very clear, positive feedback that the ampoule has been properly activated.



The packaged items, ready to be sterilised, are placed inside the supplied sterilisation bag. The cabinet holds the space around the sterilisation bag below room pressure. This ensures that any leak is from the room into the cabinet. Any EO released through the bag is taken to the outside through the exhaust tube at the back of the cabinet.

The bag volumes vary in size, giving great flexibility in the volume, size and shape of items that can be sterilised.

**Large Ampoule and Bag Refill Kit**

**AN7916**

**\$850**

consists of:

- 14 x 18cc Ampoules
- 14 x 40 Litre Sterilisation Bags
- Probe Strap

**Andersen Dosimeter Chemical Integrator**



Time, temperature and EO concentration are essential to successful sterilisation. Dosimeter integrators present visual assurance that all these parameters have been fulfilled. Dosimeters change colour in proportion to the dose of sterilant they contact. View results immediately after sterilisation. No laboratory culture is necessary.

Use one Dosimeter per cycle.

**Features:**


- Read results immediately
- Reliable
- Easy to use for every sterilisation load

**Contact Us for Reference Sites**

**AN87**

**\$125**  
Box of 25

**EO Indicator Stickers**



The EO Sterilisation Indicator Sticker is designed to provide immediate visual results as to whether or not Ethylene Oxide gases have successfully penetrated the sterilisation load.

**\$50**  
500 Stickers  
**\$100**  
1000 Stickers



### Endoscope Cleaning Brushes

Various sized brushes

- 1.0-1.5mm Instrument Channel/240cm Working Length
- 2.0-2.8mm Instrument Channel/220cm Working Length
- 2.0-2.8mm Instrument Channel/350cm Working Length

**\$10**



### Matrix Wipes

Unique biofilm removing detergent wipes for pre-cleaning endoscopes and other medical devices. Ready to use and conveniently packaged wipes. Contain MatrixT biofilm removing solution (a patented formulation of surfactants specifically designed for removing biological soils and biofilms). Apetured wipes to ensure sufficient density for removing organic soiling.

**\$15**



### Matrix Biofilm Remover Cleaning Solution

Matrix is a breakthrough in validated cleaning protocols. More effective than enzymatic cleaners, Matrix is the only proven solution designed to remove bacterial biofilms in which both bacteria and other microorganisms may be imbedded on flexible and rigid endoscopes and accessories which have accumulated over time due to insufficient cleaning.

**\$200**

5 Litre

**Impacts 100% on Viable Bacteria to Provide Absolute Assurance That All Your Endoscopes and Instruments are Optimally Cleaned.**



### Medizyme Multi-Enzyme Cleaning Solution

Medizyme **Multi-Enzyme Cleaner** is totally compatible with all flexible and rigid endoscopes due to its neutral PH. Provides safe and effective cleaning against proteins, carbohydrates, lipids and mucopolysaccharides.

**\$25**

500ml

Also provides a combination of better detergency and good rinsing ability to ensure the best possible results.

**\$150**

5 Litre



### Endozyme Sponge

Endozyme Sponges are uniquely contoured to fit around the outside sheath of scopes and are pre-saturated with Endozyme Bio-Clean to provide a quicker, more thorough pre-cleaning.

**\$10**

Each



### OPAL Disinfectant

OPAL™ is an Ortho-Phthalaldehyde (OPA) based **high level liquid disinfectant** that features a reuse period of up to 21 days and a fast turnaround of 6-minute contact time at 20°C. Ready-to-use, no mixing or activation required. Effective against a wide array of micro-organisms. Suitable for use on all flexible and rigid endoscopes and accessories.

**\$150**

5 Litre



### OPAL Test Strips

Opal® Test Strips are chemical indicators for use in determining whether the concentration of OPA, the active ingredient in OPAL®, is above or below the Minimum Recommended Concentration of 0.35%. It is recommended the OPAL® solution be tested before each disinfection cycle.

**\$60**

50 Test Strips

- Easy to use - 3 second dip time
- Results in 90 seconds
- Colour reference chart to confirm concentration amount. Large, easy to read indicator pad



### Mr Clear Anti-Fog Solution


Convenient 6ml Anti-Fog solution which eliminates condensation at the distal end of rigid endoscopes/telescopes. Non-abrasive sponge allows for effective cleaning without distal lens damage.

**\$15**



**Austvet Endoscopy is committed to:**

- Providing the highest level of Endoscopic service and repair
- Providing the fastest turnaround time
- Maximising our customers' satisfaction and productivity
- Providing training and technical support
- **NEW General Instrument Repair Service**  
including tightening, realigning, wire brushing, sharpening, polishing, repinning, insulation testing, soldering and modifications to general and endoscopic instruments



**In-House Service & Repair.  
Don't Take Risks With Your Precious Endoscopy Equipment!**





AUSTVET ENDOSCOPY



**Steriflex**  
SAFE STERILISATION SOLUTION

ARTG  
Listed



## Effective & Efficient EO Sterilisation for the 21<sup>st</sup> Century

Austvet Endoscopy in conjunction with Steriflex are dedicated to the safe and effective sterilisation of delicate medical/veterinary devices. We are the sole Australian and New Zealand Distributors for the Andersen EO GAS Flexible Chamber Technology.

EO Gas not only offers the highest level of sterility but does so with negligible wear and tear to the instruments.

This not only saves money and time but provides you with peace of mind knowing that your equipment is sterile and that you are also helping the environment by reusing instruments and saving on landfill and waste.

EO (Ethylene Oxide) Gas is the only sterilisation modality that can cater for the sterilisation of the widest range of Endoscopes available today.

**Refer to pages 32 & 33.**

We are also proud to represent

