



**AUSTVET ENDOSCOPY**

**PRODUCT CATALOGUE**

**2021**

**Est. 2004**



## **The Future of Veterinary Endoscopy**

Please don't hesitate to call any time with any questions or queries  
or come and visit our showroom soon:

---

**Unit 5 / 14 Lionel Rd, Mt Waverley VIC 3149**

**Maria: 0409 383 327 | Ennio: 0408 053 100 | Wally: 0438 084 418**

**[sales@austvetendoscopy.com.au](mailto:sales@austvetendoscopy.com.au)**

**[www.austvetendoscopy.com.au](http://www.austvetendoscopy.com.au)**



Distributorships	3
Olympus HD Video Endoscopy Systems	4
Portable HD Video Endoscopy Systems	8
Retrieval Devices	10
Therapeutic Devices	11
Fluid Administration & Therapy	12
Endoscopic Accessories	13
Light Sources	14
Endoscopic Video Cameras / Insufflators & Accessories	15
Endoscopic Surgical Tower	16
Arthroscopic Fluid Management	17
Otoscopy/Rhinology	18
Portable HD Video Rigid Endoscopy System	19
Arthroscopy	20
Laparoscopy	26
Caiman Vessel Sealing Device	29
Trocars & Cannulas	30
ETO Gas Sterilisation	34
Service & Repair	36
Cleaning Solutions	38
Microscope	39

**DISCLAIMER:** Please note that prices and availability are correct at the time of publishing but may be subject to change without notice.

**\*All prices exclude GST and freight**



**Austvet Endoscopy is the exclusive distributor of Olympus Endoscopy products in Australia and New Zealand.**

**Recognised as the world's leading supplier of medical endoscopy equipment, we are proud to offer to the veterinary industry the Olympus range of products.**



Olympus Australia represents many years of planning and sets new standards in delivering a range of quality services and products to the Australian market. Olympus Australia was formed in 1997 and has a subsidiary company Olympus New Zealand, formed in 1998 .

Today, Olympus sells and services endoscopy and surgical equipment to the medical industry, microscopes used for research, medical and industrial purposes and remote visual testing equipment used by industry.

The medical instrument division encompasses a broad range of endoscopy and surgical instruments, with particular emphasis in the fields of:

Gastroscopy, colonoscopy, bronchoscopy and the entire spectrum of endoscopic surgery, including laparoscopy, cystoscopy and arthroscopy.



**We are also proud to represent Andersen, Aohua, B. Braun & Dr. Fritz.**

**This allows Austvet Endoscopy to provide the complete endoscopy solution to our customers.**





## OLYMPUS®

The Olympus 190 Video System is cutting-edge technology for unsurpassable visualisation.

### Endoscopy Video System

#### Includes:

- Olympus CV-190 Video Processor with Keyboard and Scope Cable (MAJ-1430)
- Olympus CLV-190 Xenon Light Source with Water Bottle
- Endoscopy Trolley with Olympus OEV-261H High Definition 26" Monitor, Fully Adjustable Monitor Arm and Endoscope Hanger

#### Features:

- Full HD imaging capabilities when used in conjunction with H-Series video endoscopes.
- Narrow Band Imaging to enhance the visibility and characterisation of healthy and abnormal tissue and other structures on the mucosal surface
- Compatible with all Olympus Endoscopes & Surgical Camera Heads

**\$27,900**

Image capture system is optional  
(see following page)







**New Endoscopy Trolley**

**\$1,600**


Trolley designed especially for Video Endoscopy.

**Features:**

- Heavy duty steel construction
- 4 tiers
- Medical grade castors with 100mm wheels
- Front castors with brakes
- Conceals all electrical wiring
- Power board 6-way, fully adjustable monitor arm and endoscope hanger

Dimensions: 550mm W x 550mm D x 1350mm H  
(Height excludes monitor)

**Monitor Optional Extra**



**New Microsoft Surface Laptop 3**

**\$3,150**

Microsoft Surface Laptop 3 with custom **HD** Image Capture Software.

**Features:**

- 13.5" display
- 2256 x 1504 screen resolution
- **HD** image and **HD** video capture capabilities
- Ability to input patient details (up to 80 characters)
- Supplied with SDI or DVI to USB converter

*Also note: Due to computer technology changing rapidly specifications and price may differ. Please enquire.*



## Specifications &amp; Compatibility Chart

Endoscope Model	Working Length	Diameter	Biopsy Channel	Angulation
URF-V	670mm	3.2mm	1.2mm	2 Way
CYF-V2	380mm	5.4mm	2.2mm	2 Way
CYF-VH	380mm	5.5mm	2.2mm	2 Way
BF-P180	600mm	4.9mm	2mm	2 Way
BF-Q180	600mm	5.2mm	2mm	2 Way
BF-P190	600mm	4.2mm	2mm	2 Way
BF-Q190	600mm	4.8mm	2mm	2 Way
BF-H190	600mm	5.5mm	2mm	2 Way
GIF-XP150N	1100mm	5.5mm	2mm	4 Way
GIF-Q160	1030mm	9.5mm	2.8mm	4 Way
GIF-N180	1100mm	4.9mm	2mm	2 Way
GIF-XP180N	1100mm	5.5mm	2mm	4 Way
GIF-Q180	1030mm	8.8mm	2.8mm	4 Way
GIF-H180	1030mm	9.8mm	2.8mm	4 Way
GIF-H190	1030mm	9.2mm	2.8mm	4 Way
PCF-H180AL	1680mm	11.8mm	3.2mm	4 Way
PCF-H190AL	1680mm	11.7mm	3.2mm	4 Way

*\*All Second Hand Endoscopes Pending Availability*



Image	Air / Water Facility	Processor CV-190	Processor CV-180	Price Second Hand
	No	✓	✓	\$9,500
	No	✓	✓	\$3,500
	No	✓	✓	\$6,900
	No	✓	✓	\$8,900
	No	✓	✓	\$8,900
	No	✓	✗	\$12,900
	No	✓	✗	\$12,900
	No	✓	✗	\$13,500
	Yes	✓	✓	\$9,500
	Yes	✓	✓	\$4,900
	Yes	✓	✓	\$12,900
	Yes	✓	✓	\$14,900
	Yes	✓	✓	\$7,900
	Yes	✓	✓	\$10,900
	Yes	✓	✗	\$14,900
	Yes	✓	✓	\$10,900
	Yes	✓	✓	\$14,900

OLYMPUS FLEXIBLE VIDEO ENDOSCOPES





- Dimensions: 48cm x 42cm x 29cm
- Weight: 17kg

Portable battery powered **HD** video-endoscopy system designed for ease of use for Veterinary surgeons performing upper airway and gastroscopic examinations.

### Endoscopy Video System

**\$15,000**

#### Includes:

- Battery Video Processor with battery life indicator (over 3 hours of battery life)
- Microsoft Surface Book with 13.5" Screen
- Custom Made Carry Case
- Easy to use **HD** Image/Video Capture Software

### VET-9215 HD Video Endoscope

**\$14,800**

#### Specifications:

- Working Length: 1500mm
- Diameter: 9.2mm
- Biopsy Channel: 2.8mm
- Air/Water Facility
- 4 Way Angulation
- Includes Leak Tester, Water Bottle and Cleaning Accessories

### VET-6010 HD Video Endoscope

**\$16,390**

#### Specifications:

- Working Length: 1000mm
- Diameter: 6mm
- Biopsy Channel: 2mm
- Air/Water Facility
- 4 Way Angulation
- Includes Leak Tester, Water Bottle and Cleaning Accessories

### VET-9830 HD Video Endoscope

**\$15,350**

#### Specifications:

- Working Length: 3000mm
- Diameter: 9.8mm
- Biopsy Channel: 2.8mm
- Air/Water Facility
- 4 Way Angulation
- Includes Leak Tester, Water Bottle and Cleaning Accessories







**AOHUA MVE-9215  
Portable Endoscope**

**When portability is A MUST**

- Working Length: 150cm
- Diameter: 9.2mm
- Biopsy Channel: 2.8mm
- 4 Way Angulation
- LED Illumination

Scope with Small Monitor

**\$9,570**



**Clax Smart Fold Cart**

**\$450**

- Can be collapsed quickly and easily at the push of a button
- Has two large platforms, supplied standard with one collapsible storage box
- Has many useful features such as foot brake, holding device for boxes, upward folding platform and removable wheels for space-saving storage
- Made of injection moulded plastic and aluminium for long lasting use

**Specifications:**




60kg Total capacity	7kg Unit weight
20kg Top load capacity	40kg Bottom load capacity
Folding Handle type	1 crate Included
890 x 550 x 1030mm Dimensions assembled	670 x 470 x 110mm Dimensions folded

Great for transporting the Portable Equine Endoscopy System



**New Oesophageal Rigid Foreign Body Removal Forceps**

Teeth face backwards to maximise grip and minimise damage to the oesophagus.

	Diameter at Jaws in mm	Length at Jaws in mm	Shaft Diameter in mm	Shaft Length in mm	Price
	13 x 9	34	7.5 x 6	640	<b>\$450</b>
	11.5 x 6	31	6.5 x 5	500	<b>\$450</b>
	5.6 x 4.2	15.5	4 x 3	400	<b>\$450</b>



“Standard” Size Most Popular

**New Flexible Retrieval Devices****Retrieval Grasping Basket**

For the extraction of round, slippery objects.

**\$330**

Paediatric: For 2.0mm > Instrument Channels, 180cm Length Basket Opening 32mm  
Adult: For 2.8mm > Instrument Channels, 230cm Length Basket Opening 30mm

**Retrieval Tripod Prongs**

For the collection of irregular shaped objects.

**\$160**

Paediatric: For 2.0mm > Instrument Channels, 180cm Length  
Adult: For 2.8mm > Instrument Channels, 230cm Length

**SnareMaster Oval Loop**

Designed with optimal stiffness to precisely capture foreign bodies. Handle shaft is marked at 5mm increments to indicate the snare throw opening.

**\$60**

Adult: For 2.8mm > Instrument Channels, 230cm Length Snare Opening 30mm  
For 2.8mm > Instrument Channels, 350cm Length Snare Opening 25mm

**SnareMaster Crescent Shape**

Designed with optimal stiffness to precisely capture foreign bodies. Handle shaft is marked at 5mm increments to indicate the snare throw opening.

**\$60**

Paediatric: For 2.0mm > Instrument Channels, 230cm Length Snare Opening 25mm

**Spider-Net Retrieval Device**

For the collection of round, blunt or otherwise hard to grasp foreign bodies. The “spider” net delivers a flexible, durable pouch eliminating the risk of accidental loss during retrieval.

**\$160**

Adult: For 2.8mm > Instrument Channels, 230cm Length Net Opening 3x6cm

**Ratigator Reusable Grasping Forceps**

Combines a rat tooth with alligator jaws into one instrument. 15mm jaw opening. Reusable & autoclavable.

**\$550**

Paediatric: For 2.0mm > Instrument Channels, 210cm Length  
Adult: For 2.8mm > Instrument Channels, 230cm Length  
For 2.8mm > Instrument Channels, 350cm Length



**New Flexible Biopsy Forceps**

	Instrument Channel	Working Length	Price
<b>Biopsy Forceps/Rat Tooth (reusable)</b>	1.2mm	115cm	<b>\$950</b>
<b>Serrated Cup Paediatric</b>	2.0mm	160cm	<b>\$50</b>
<b>Oval Cup Paediatric</b>	2.0mm	160cm	<b>\$20</b>
<b>Serrated Cup w Spike Paediatric</b>	2.0mm	160cm	<b>\$50</b>
<b>Serrated Cup Adult</b>	2.8mm	230cm	<b>\$20</b>
<b>Serrated Cup w Spike Adult</b>	2.8mm	230cm	<b>\$20</b>
<b>Oval Cup w Spike Adult</b>	2.8mm	230cm	<b>\$20</b>
<b>Oval Cup Adult</b>	2.8mm	230cm	<b>\$20</b>
<b>Oval Cup Adult</b>	2.8mm	400cm	<b>\$120</b>



**Flexible Cytology Brushes**

Designed to enhance the efficiency of bronchoscopic and gastrointestinal examinations and improve collection of cytology samples. Features optimal brush hardness for better sample collection and thinner bristles for improved cell collection.

	Instrument Channel	Working Length	Brush Size	Price
<b>Bronchial</b>	1.2mm	115cm	2mm	<b>\$50</b>
	2.0mm	120cm	3mm	<b>\$30</b>
	2.0mm	120cm	2mm	<b>\$30</b>
<b>Gastrointestinal</b>	2.0mm	200cm	3mm	<b>\$30</b>
	2.8mm	200cm	3mm	<b>\$30</b>



**Injector Needle**

Designed to provide proper flexibility and rigidity to ensure safe piercing and excellent puncture performance. Ideal for injecting Formalin to treat Ethmoid Haematomas in horses. Assorted needle lengths and gauges. Large range of needle gauges and lengths for paediatric and adult scopes, please enquire

Instrument Channel	Working Length	Needle Diameter	Needle Length	Price
2.0mm	200cm	23G	5mm	<b>\$70</b>
2.0mm	165cm	21G	6mm	<b>\$70</b>
2.8mm	230cm	23G	6mm	<b>\$70</b>
2.8mm	240cm	21G	8mm	<b>\$70</b>



**Olympus EZDilate Inflation Device**

**\$150**

Used while performing balloon dilation procedures to inflate the balloon, monitor the pressure within the balloon and deflate the balloon. Easy to read luminescent dial and quick release lock lever mechanism. Maximum pressure of 30 atmospheres.



**Olympus EZDilate Multi-Stage Endoscopic Balloon Dilators**


**\$350**

Designed for oesophageal dilation under direct endoscopic vision. Low compliance balloon material produces dependable, controllable radial dilation. Low profile balloon design assures smooth accessibility into difficult strictures.

Assorted Sizes: For 2.8mm > Instrument Channels, from 6mm to 20mm Diameter

<b>Inflation diameter (mm)</b>	6, 7, 8	8.5, 9.5, 10.5	11, 12, 13	13.5, 14.5, 15.5	16, 17, 18	18, 19, 20
<b>Inflation pressure (atm)</b>	2, 5.5, 9.5	2, 5.5, 9.5	2, 3.5, 5.5	2.5, 4, 6	2.5, 4, 5.5	2.5, 3.5, 5
<b>Balloon length (mm)</b>	65	65	55	55	55	55





### Mila Double/Triple Lumen Flushing Catheter

Flushing catheter for guttural pouch lavage. Can also be used as a sampling catheter for sterile specimen retrieval via an endoscope. The set contains a glycol plug in the outer catheter and an inner catheter for retrieval of the sterile specimen.

**1.8mm Diameter/250cm Length/Double Lumen**

**\$80**

**2.3mm Diameter/190cm Length/Triple Lumen**

**\$80**



### Equine Uterine Flushing Catheter

Used for flushing and irrigation, administration of medications, embryo collection and transfer. Autoclavable. Sterile.

**12mm Diameter/80cm or 150cm Length**

**\$132**



### Endoscopic BAL Small Animal/Equine Single Lumen Flushing Catheter

An endoscopically guided catheter for lung lavage or introduction into the guttural pouch to flush out pouch debris. Visual indicators at 2.5cm, 5cm and 7.5cm from the distal tip to aid in introduction. Sterile.

**1.8mm Diameter/120cm Length**

**\$30**

**1.8mm Diameter/250cm Length**

**\$50**

**2.3mm Diameter/220cm Length**

**\$15**

**2.3mm Diameter/400cm Length**

**\$20**



### Silicone Stomach Tube

Three strategically placed fenestrations at the distal tip provide improved safety and effectiveness. Includes syringe and funnel adaptor. Autoclavable. Non-Sterile. (Internal dimensions listed).

**4.5mm/91cm Length or 5.0mm/91cm Length**

**\$80**

**5.0mm/183cm Length or 6.5mm/150cm Length**

**\$80**

**6.5mm/300cm Length or 9.5mm/300cm Length**

**\$115**

**11mm/300cm Length or 16mm/305cm Length**

**\$150**

**10mm/300cm Length (Standard Line)**

**\$70**



### Broncho-Alveolar Lavage Catheter (Silicone) BAL TUBE

The BAL catheter is used for lung lavage of selected segments of the lung fields in horses. Equipped with a distal balloon that can be inflated to prevent retrograde fluid flow and a proximal two-way fitting allows simultaneous connection for both irrigation and aspiration. Sterile and Autoclavable

(The plastic proximal two-way fitting and three-way stopcock are not autoclavable).

**10mm Diameter/240cm Length**

**\$170**

**10mm Diameter/300cm Length**

**\$180**

**Three-Way Replacement Stopcock**

**\$20**





**Standard Single Scope Hanger**

Requires only two screws for mounting.  
Dimensions: 145mm W x 135mm D

**\$150**



**Olympus, Aohua and Pentax Water Bottles**

Pentax Water Bottle to suit all Pentax Endoscopes  
Olympus Water Bottle Autoclavable for 140/160/180 Series  
Olympus Endoscopes  
Olympus Water Bottle Autoclavable for 100 Series  
Olympus Endoscopes

**\$500**

**\$280**

**\$190**

**Second Hand Bottles Also Sold Pending Availability**



**Olympus, Aohua and Pentax Leakage Testers**

Designed to detect and identify leaks from your flexible endoscopes. Stops potential fluid invasion thereby preventing the possibility of costly repairs, cross infection and downtime.

**\$300**

**ALL flexible endoscopes must be leak tested prior to EVERY submersion to ensure YOU, YOUR PATIENTS AND YOUR SCOPES are protected.**



**Askir Endoscopic Suction Pump & Accessories**

Compact, high performance unit designed for both endoscopic aspiration and surgical suction. Unit is supplied complete with polycarbonate autoclavable 1 litre jar with overflow valve and antibacterial filter. Pump is equipped with aspiration regulator and vacuum indicator. Designed for almost continuous use up to 120 minutes. Maximum suction flow 30 litres per minute. Portable and lightweight.

**\$550**

**Accessories:**



Busse Mucus Specimen Trap  
80cc

**\$10**



Suction Handpiece Frazier  
2.6mm, 3.3mm or 4.0mm

**\$25**



Suction Handpiece Yankauer 29cm  
Standard or Paediatric

**\$35**







**NEW AE-L100 Endoscopic LED Light Source (No Air Pump)**

**\$1,200**

Ideal for use with all rigid endoscopes including Arthroscopes, Laparoscopes, Otosscopes, Rhinoscopes and Cystoscopes.

**Features:**

- Extremely compact (30cm x 20cm x 10cm)
- Expected lifetime of LED is over 40,000 hours
- Extremely bright light
- Supplied with light guide cable



For Rigid Endoscopes

**LED Portable High-Intensity Light Handles for Rigid and Flexible Endoscopes Now with Variable Brightness**

LED Portable Light Handles are an easy option for quick examinations in the consult room or for fieldwork. They feature high-intensity illumination and are a great light source backup. Handles attach to most rigid and flexible endoscopes.



For Flexible Endoscopes

**Features:**

- Lithium-Ion rechargeable battery
- Extremely compact
- High-intensity white light with variable brightness

**Kit Includes:**

- LED light source handle
- Rechargeable Lithium-Ion battery
- Battery charger

LED Light Handle for Rigid Endoscopes **\$330**



LED Light Handle for Flexible Endoscopes **\$330**




**Olympus Endo LED Portable Light Handle for Rigid Endoscopes**

**\$1,350**



- Colour temperature like Xenon light
- High brightness over entire LED lifetime of 100,000 hrs
- Convenient battery concept using standard AA or rechargeable batteries and a battery status indicator
- Ergonomic design for rigid endoscopes for more comfort
- Fully immersible, washer and disinfectant compatible
- 2-3 Hour Duration



**Battery Indicator**

For more safety, the battery status indicator shows green for fully charged, yellow for partially discharged or flashing yellow for fully discharged



**Stryker 1288 HD**

**\$3,990**

Late model Stryker Camera with HD imaging capabilities. Camera head clips directly onto all rigid endoscopes to provide visualisation of the surgical site across all endoscopic applications. In addition, with a video adaptor, the camera can be connected to all Olympus and Pentax flexible endoscopes.



**New Sterile Camera Sleeve with Telescope Adaptor**

**\$30** Each  
or  
**\$1,080**  
Box of 40

Disposable surgical sterile video camera sleeve/cover with enclosed adaptor for telescope connection. Safe, quick and easy application. Latex free. Offers 100% sterility without chemicals to reduce risk of infection. Excellent equipment protection.

Dimensions: 14cm W x 285cm L



**Olympus UHI-3 35 Litre Insufflator**

**\$690**

Late model Olympus insufflator which features preset pressure and flow rate, leakage warning display and overpressure warning display and alarm. Supplied with hose for gas bottle. Designed to support a wide range of Laparoscopic insufflation applications.



**New Laparoscopic Insufflator Tubing**

**\$10**

Disposable CO2 Insufflator Filter with 3m Tubing. Has universal fitting however any additional connectors required can be supplied at no charge. Filtration area 24cm<sup>2</sup> for added protection and greater gas flow. Sterile and Latex free.



**New Verres Needles**

Verres Needles equipped with a spring loaded obturator to insufflate the abdomen in Laparoscopic surgery.

100mm	<b>\$145</b>
150mm	<b>\$163</b>
200mm	<b>\$210</b>





**Stryker 1288 Surgical Endoscopic Tower  
Ideal for ENT, Arthroscopy, Laparoscopy and  
Otoscopy.**

**System Includes:**

- Stryker 1288 HD Video Camera Console and Camera Head
- Stryker L9000 LED Light Source (Without Air Pump)
- Stryker 40L Hi-Flow Insufflator
- Stryker Vision Elect HD 26" Flat Panel Monitor with Articulating Arm
- Stryker Light Guide Cable
- Stryker Endoscopy Trolley with Drawer and Glass Doors

**\$8,950**





**New Linvatec 10k Pump. The Simplest High-Performance Pump in Arthroscopy.**

**Ease of Use**

- Intuitive one-step/one-hand cassette loading
- Quick Setup and priming
- Integrated pole mount works with any standard I.V. pole
- Compatible with video cart shelf or any flat surface

**Superior performance**

- Immediate response
- High flow
- Maintains distension under even the most demanding situations
- Robust autoclavable remote control

**Specifications**

- Maximum set pressure 150mmHg
- Maximum flow rating 2000ml/min
- Weight: 3kg
- Dimensions: 16.5cm tall, 20.3cm deep and 19.7cm wide
- Voltage 90–240v



**\$4,600**

**10k100 Endoscopy Tube Set**

**\$650**  
box of 10



**New Dr Fritz Pressure Infusion Cuff**

500ml - 1 Litre (with manometer)	<b>AR-08-100N</b>	<b>\$390</b>
3 - 5 Litre (with manometer)	<b>AR-08-300</b>	<b>\$690</b>



**New Standard Veterinary Rhinoscope Set****\$5,190**

Rhinoscope set supplied with two sheaths, with the larger sheath allowing both biopsy and irrigation.

Features universal eyepiece for compatibility with all Endoscopic video cameras.

Multi-purpose rigid endoscope which can also be used for:

**Arthroscopy**

**Avian Endoscopy**

**Cystoscopy**

**Laparoscopy**

**Thoracoscopy**

**Vaginoscopy**

**System Includes:**

- Rhinoscope 2.7mm 30 Deg (190mm Working Length)
- 3.7mm/11 Fr Sheath (no biopsy facility – only irrigation)
- 4.7mm/14 Fr Sheath (biopsy & irrigation)
- High-intensity LED Portable Light Handle (with variable brightness)
- Cleaning Brush & Custom Carry Case
- Supplied with your choice of 1.6mm Biopsy or Grasping Forceps (see page 21)

Otoscope not included (please refer to below)

**New Standard Veterinary Otoscope Set****\$1,840**

Complete, all-in-one veterinary otoscope. Ideal for quick examinations in the consult room.

Large field of view and depth of field.

Also used for foreign body removal and/or biopsy and quick diagnostics.

**Set Includes:**

- Otoscope 5.5mm 0 Deg (100mm Working Length)
- High-intensity LED Portable Light Handle (with variable brightness)
- Two Way Stopcock (Removable)
- Cleaning Brush & Carry Case (refer to image of carry case on page 19)
- Supplied with your choice of 1.6mm Biopsy or Grasping Forceps (see page 21)



**Essential Instrument for Every Clinic**



**New All-in-One Portable Video System (AE-1080 HD)**

**\$7,100**

This all-in-one portable **HD** Video System is a complete, user friendly solution for all types of rigid endoscopy for the practice as well as the mobile vet.

**Features:**


- 15.5" full **HD** display
- 3 metre long camera and light lead cables
- Image/video capture facilitated from camera head
- Multiple video outputs for slave monitors

**System Includes:**

- 15.5" full **HD** display
- Video processor
- Camera head
- LED light source
- Light Lead
- Custom carry case



**48404**



Ultra slim 1.6mm/5fr semi flexible biopsy forcep. Ideal for biopsy of ear, nose or throat.

**New 1.6mm/5fr Semi Flexible Biopsy/Retrieval Forceps**


**\$300**

Will suit rigid ENT/Otosopes with instrument channels 1.8mm and above.

Working length: 270mm.




**48405**



Ultra slim 1.6mm/5fr semi flexible grasping forcep. Ideal for removal of foreign bodies from ear canals and nasal cavities.



**48405-B**




Ultra slim 1.6mm/5fr semi flexible dual purpose biopsy/retrieval forcep.





2.7mm

**A70963A** Telescope,  
**\$3,043** Autoclavable,  
 30° direction of view,  
 Wide angle field of view,  
 70 mm working length



**WA70992A** Trocar tube,  
**\$1,060** Rotatable,  
 1 Stopcock,  
 Standard flow,  
 for 70 mm telescopes,



**A70978A** Trocar spike,  
**\$262** triangular,  
 for 70 mm telescopes



**A70979A** Trocar spike,  
**\$262** blunt,  
 for 70 mm telescopes




**A00507A** Sealing, set  
**\$59** Set of spare sealing rings:  
 Contains 2 large sealings  
 (Transparent) and 5 small sealings (black).



**A0460** Light-guide adapter,  
**\$58** for OES 4000 light-guide cables

**00332** Light-guide adapter,  
**\$23** for Wolf light-guide cables


**A5333** Eyepiece cup,  
**\$96** for telescopes



4mm



**A70941A** Telescope,  
**\$3,167** Autoclavable,  
 30° direction of view,  
 115° field of view,  
 160mm working length



**WA70952A** Trocar tube,  
**\$1,060** Rotatable,  
 1 Stopcock,  
 Standard flow




**A70954A** Trocar spike,  
**\$270** Triangular



**WA70953A** Trocar tube,  
**\$1,060** Rotatable,  
 2 Stopcocks,  
 Standard flow



**A70955A** Trocar spike,  
**\$270** Blunt




**New Standard Line 3mm ENT/Arthroscope**

Arthroscope supplied with 4mm Cannula and Sharp and Blunt Trocar.

**\$1,500**

**Specifications:**

- Working Length: 175mm (with cannula 135mm)
- Diameter: 3mm
- Direction of View: 30 Degree



**New Standard Line 4mm ENT/Arthroscope**

Arthroscope supplied with 5mm Cannula and Sharp and Blunt Trocar.

**\$1,500**

**Specifications:**

- Working Length: 175mm (with cannula 135mm)
- Diameter: 4mm
- Direction of View: 30 Degree





### New Dr Fritz Small Animal Arthroscopy Hand Instruments

Instruments provide easy handling and an ergonomic design with handles. All instruments are specifically designed for Small Animal Arthroscopy.



#### AR-17-03S

Bayonet Knife with Handle, 2.7mm

\$350



#### AR-17-05S

Banana Knife with Handle, 2.7mm

\$350



#### AR-17-12S

Mikro Picking Knife with Handle, 2.7mm

\$350



#### AR-17-19S

Hook Knife with Handle, 2.7mm

\$350



#### AR-17-06S

Curette Curved with Handle, 2.7mm

\$350



#### AR-17-07S

Ring Curette with Handle, 2.7mm

\$350



#### AR-17-11S

Meniscutome Smilie with Handle, 2.7mm

\$350



#### AR-17-09S

Hook Probe with Handle, 2.2mm

\$350



#### AR-17-13

Special Elevator with Handle, 2.7mm (Metal)

\$350



#### AR-17-22

Coronoid Chisel with Handle, 3/4mm Working Length 20cm (Metal)

\$350



#### I-17-822 (SET)

Milling Drill  
2.2mm, Working Length 12cm,  
Rotating Handpiece and Long Lasting  
Drill Steel (Dismantable)

\$500



Replacement Drill Steel (for Milling Drill)

\$165



#### S110-1935

Protection & Transportation Sheath for  
Small Arthroscopes 12cm, Autoclavable

\$40

**New Dr Fritz Small Animal Arthroscopy  
Punches, Forceps & Scissors**



**AS-17-521** **\$1,400**  
Mini Scissors for Arthroscopy, 2mm,  
Working Length 11cm



**AT-19-231** **\$1,400**  
Mini Rongeur & Biopsy Forceps, 2.2/2.7mm,  
Working Length 11cm



**AS-17-632** **\$1,400**  
Universal Rongeur & Grasping Forceps, 2.7mm,  
Working Length 12cm



**AS-17-662** **\$1,400**  
Alligator Grasping Forceps with Ratchet, 2.7mm,  
Working Length 12cm



**AS-17-621** **\$1,400**  
Arthroscopic Hook Scissor, 2.7mm,  
Working Length 12cm



**AS-17-672** **\$1,400**  
Arthroscopic Punch, 2.7mm,  
Working Length 12cm



**AS-17-641** **\$1,400**  
Biopsy & Spoon Forceps, 2.7mm,  
Working Length 12cm

**New Dr Fritz Small Animal Working Cannula Set**

This system allows the specific exchange of different working cannulas. Beginning with the smallest diameter, the vet can always change to a bigger diameter by using the changing rods supplied. Cannulas are specifically designed for Small Animal Arthroscopy.



**I-17-755** **\$890**  
**Working Cannula Set**

Consists of:

- 3 Cannulas for Elbow 2.3mm, 2.9mm and 3.5mm, Working Length 3cm – 5.5cm
- 2 Cannulas for Shoulder 2.9mm and 3.5mm, Working Length 7cm – 8cm
- 2 Sharp Trocars
- 2 Changing Rods, 2.3mm and 2.8mm





### New Dr Fritz Equine Arthroscopy Hand Instruments

Instruments have been specifically designed for Equine Arthroscopy.

	<b>I-18-100</b> Osteotome Curved with Handle, 4mm	<b>\$220</b>
	<b>I-18-110</b> Osteotome Straight with Handle, 4mm	<b>\$220</b>
	<b>I-18-520</b> Hook Probe with Handle, 3mm, Working Length 12cm	<b>\$220</b>
	<b>I-18-521</b> Special Knife with Handle, 3mm, Break Safe Special Cut	<b>\$260</b>
	<b>I-18-963</b> Hook Scissors, Straight, Single Action Jaws, 4.5mm, Working Length 18cm	<b>\$760</b>
	<b>I-18-601</b> Curette Straight with Handle, Oval Cup 6 x 4mm	<b>\$220</b>
	<b>I-18-621</b> Curette, 35 Degree Oblique Left, Oval Cup, 6 x 4mm	<b>\$220</b>
	<b>I-18-622</b> Curette, 35 Degree Oblique Right, Oval Cup, 6 x 4mm	<b>\$220</b>
	<b>I-18-518</b> Egress Cannula with Sharp Trocar and Stopcock, 2.7mm, Working Length 7cm	<b>\$270</b>
	<b>I-18-519</b> Egress Cannula with Sharp Trocar and Stopcock, 4mm, Working Length 7cm	<b>\$270</b>
	<b>I-18-635</b> Special Egress Curette Straight, Oval Cup, 5.5 x 3mm and Leur Lock (f) Suction Connector, Working Length 12cm	<b>\$410</b>
	<b>I-18-536</b> Special Elevator for Coffin Joint, Working Length 8cm	<b>\$235</b>
	<b>I-18-517</b> Elevator Curved with Handle, 6mm	<b>\$275</b>
	<b>I-18-534</b> Special Elevator for Sesamoide Fractures, Curved, Long 4mm	<b>\$315</b>



**New Dr Fritz Equine Arthroscopy Line Instruments**

**“Standard Line” Rongeur**

**I-18-6611**

**\$950**

Straight, Fenestrated, Oval Cup 4 x 8mm, Working Length 11cm

**I-18-662**

**\$950**

Angled Upwards, Fenestrated, Oval Cup 4 x 8mm, Working Length 11cm



**“Power Line” Rongeur**

**I-18-973**

**\$1,190**

Straight, Heavy Durable Handle, Oval Cup 5 x 10mm, Working Length 15cm

**I-18-974**

**\$1,190**

Straight, Heavy Durable Handle, Fenestrated, Small Tip, Oval Cup 5 x 10mm, Working Length 15cm



**New Dr Fritz Equine Arthroscopy Ergonomic Instruments**

**Straight Rongeur**

**I-18-7613**

**\$1,190**

Straight, Ergonomic Handle, Fenestrated Cup 4 x 8mm, Working Length 12cm



**Hook Grasper “Chip Catcher” Rongeur**

**I-18-6414**

**\$1,190**

For Large Fragments, Ratchet, Working Length 12cm



**New Ferris Smith Equine Rongeurs**

**\$990**

Premium quality German made Rongeurs. Features straight and angled upward tip with 180mm working length. Lifetime warranty.

**Sizes Available:**

- 3mm W x 10mm L Tip (Straight & Angled Upward)
- 4mm W x 10mm L Tip (Straight & Angled Upward)
- 5mm W x 10mm L Tip (Straight & Angled Upward)
- 6mm W x 10mm L Tip (Straight only)



### New Dr Fritz Needle Holder with Handle

Needle self positioning, axial handle.

Small Animal: **L-19-562V** **\$1,090**

Left Turned 5mm,  
Working Length 22cm

**L-19-563V** **\$1,090**

Angled Upwards 5mm,  
Working Length 22cm

Equine: **L-19-591V** **\$1,090**

Straight 5mm,  
Working Length 43cm

### New Dr Fritz Small Animal 5mm Hand Instruments

Instruments are easily dismantable and rotatable through 360 Degrees. Each Instrument is supplied complete with jaw insert, sheath and ratchet handle.



**531.61H01V** **\$750**

Biopsy Forceps, Double Action Jaws, With Spikes  
5mm, Working Length 33cm



**541.61H01V** **\$750**

Biopsy Forceps, Double Action Jaws, With Spikes  
5mm, Working Length 43cm



**531.01H05V** **\$750**

Mini Metzenbaum Scissors, Curved Left, Insulated  
5mm, Working Length 33cm



**541.01H05V** **\$750**

Mini Metzenbaum Scissors, Curved Left, Insulated  
5mm, Working Length 43cm



**531.22H11V** **\$750**

Grasping Forceps Universal, Insulated  
5mm, Working Length 33cm




**541.22H11V** **\$750**

Grasping Forceps Universal, Insulated  
5mm, Working Length 43cm



**531.40H01V** **\$750**

Dissecting Forceps Maryland, Turned, Insulated  
5mm, Working Length 33cm



**541.40H01V** **\$750**

Dissecting Forceps Maryland, Turned, Insulated  
5mm, Working Length 43cm




**531.30H12V** **\$750**

Grasping Forceps Claw, Straight, Metal  
5mm, Working Length 33cm



**541.30H12V** **\$750**

Grasping Forceps Claw, Straight, Metal  
5mm, Working Length 43cm



**531.27H11V** **\$750**

Grasping Forceps Babcock  
5mm, Working Length 33cm

**New Dr Fritz Knot Tier**

**L-19-534**

**\$390**

Standard Version, 5mm, 29cm  
 Standard Version, 5mm, 40cm  
 Easy filling, suture size up to 1 metric  
 Easy handling, easy to tighten



**New Dr Fritz Injection and Puncture Cannula**

Interchangable needle (supplied 1mm x 26mm) with protection outer tube and inside tube with Luer Lock connector.

**\$330**

Small Animal: **L-33-210** 5mm, Working Length 33cm  
 Equine: **L-47-210** 5mm, Working Length 47cm  
 Spare Needle **L-210** 1mm x 26mm (2 Pack)

**\$70**



**New Dr Fritz Equine 10mm Hand Instruments**

Instruments are easily dismantable and rotatable through 360 Degrees. Each Instrument is supplied complete with jaw insert, sheath and ratchet handle.

**141.11H05V**

**\$950**

Mayo Scissors, Curved Left, Insulated  
 10mm, Working Length 45cm



**141.48H11V**

**\$950**

Dissecting Forceps Overholt, Turned, Insulated  
 10mm, Working Length 45cm



**141.26H11V**

**\$950**

Grasping Forceps Babcock,  
 10mm, Working Length 45cm



**141.25H11V**

**\$950**

Grasping Forceps Universal,  
 Double Action Jaws, Insulated  
 10mm, Working Length 45cm



**141.49H11V**

**\$950**

Grasping Forceps Atraumatic, Insulated  
 10mm, Working Length 45cm



**141.32H11V**

**\$950**

Grasping Forceps Claw, Double Action Jaws, Metal  
 10mm, Working Length 42cm






### New Olympus WA29042B AI Endoscope with Sterilising Tray

Laparoscope/Endoscope supplied with Stopcock with Two Irrigation Ports, Instrument Tray and Sealing Caps. Angled ocular eyepiece. Ideal for Reproduction, ENT and Urology. Fully autoclavable.

**\$9,460**

#### Specifications:

- Working Length: 430mm
- Diameter: 8.6/9.8 Fr
- Instrument Channel: 6.4 Fr
- Direction of View: 7 Degree



### New Standard Line 11mm Operating Laparoscope

Ideal for single Port Laparoscopic Procedures. Cannula and Trocar optional extra.

**\$4,900**

#### Specifications:

- Working Length: 312mm
- Diameter: 11mm
- Instrument Channel: 6.1mm



### Olympus/Stryker 5mm Laparoscope

Ideal for Endoscopic Artificial Insemination in sheep and Small Animal/Equine Laparoscopic applications. Cannula and Trocar optional extra.

**\$2,000**

#### Specifications:

- Working Length: 300mm
- Diameter: 5mm
- Direction of View: 0 Degree or 30 Degree Available



### Olympus/Stryker 10mm Laparoscope

Ideal for Endoscopic Artificial Insemination in sheep and Equine Laparoscopic applications and equine dentistry. Cannula and Trocar optional extra.

**\$1,600**

#### Specifications:

- Working Length: 300mm
- Diameter: 10mm
- Direction of View: 0 Degree or 30 Degree Available



### Olympus/Stryker 10mm Extra Long Laparoscope

Ideal for Equine Laparoscopy or Dentistry. Cannula and Trocar optional extra.

**\$3,000**

#### Specifications:

- Working Length: 450mm
- Diameter: 10mm
- Direction of View: 0 Degree or 30 Degree Available



# AESCULAP® CAIMAN®

ADVANCED BIPOLAR SEAL AND CUT TECHNOLOGY

A multipurpose vessel sealing instrument for Laparoscopic and Open surgery.

- Unique, tip first closure with patented floating hinge for vessel encapsulation
- Seals vessels up to 7mm Ø
- 26.5mm seal length, 23.5mm cut length
- 5mm Ø (Lap) and 12mm Ø (Open) in 24cm, 36cm and 44cm (articulating)

**New, compact Seal and Cut generator with warranty and a readily available supply of hand pieces**



**Caiman GN200 Generator \$7,600**



PL738SU	CAIMAN SEALING AND CUTTING INSTRUMENT	5mm	24cm	NON ARTICULATED	<b>\$580</b>
PL740SU	CAIMAN SEALING AND CUTTING INSTRUMENT	5mm	36cm	NON ARTICULATED	<b>\$580</b>
PL742SU	CAIMAN SEALING AND CUTTING INSTRUMENT	5mm	44cm	NON ARTICULATED	<b>\$580</b>
PL739SU	CAIMAN SEALING AND CUTTING INSTRUMENT	5mm	24cm	ARTICULATED	<b>\$660</b>
PL741SU	CAIMAN SEALING AND CUTTING INSTRUMENT	5mm	36cm	ARTICULATED	<b>\$660</b>
PL743SU	CAIMAN SEALING AND CUTTING INSTRUMENT	5mm	44cm	ARTICULATED	<b>\$660</b>
PL730SU	CAIMAN SEALING AND CUTTING INSTRUMENT	12mm	24cm	ARTICULATED	<b>\$720</b>
PL731SU	CAIMAN SEALING AND CUTTING INSTRUMENT	12mm	44cm	ARTICULATED	<b>\$720</b>
PL774SU	CAIMAN MARYLAND SEALING AND CUTTING INSTRUMENT	5mm	12.5cm	NON ARTICULATED	<b>\$577</b>
PL775SU	CAIMAN MARYLAND SEALING AND CUTTING INSTRUMENT	5mm	17cm	NON ARTICULATED	<b>\$577</b>
PL770SU	CAIMAN MARYLAND SEALING AND CUTTING INSTRUMENT	5mm	36cm	NON ARTICULATED	<b>\$577</b>
PL772SU	CAIMAN MARYLAND SEALING AND CUTTING INSTRUMENT	5mm	44cm	NON ARTICULATED	<b>\$577</b>
PL771SU	CAIMAN MARYLAND SEALING AND CUTTING INSTRUMENT	5mm	36cm	ARTICULATED	<b>\$1,025</b>
PL773SU	CAIMAN MARYLAND SEALING AND CUTTING INSTRUMENT	5mm	44cm	ARTICULATED	<b>\$1,025</b>





Ø 5mm

**Sealing Units**

Reusable



**EK080P \$350**

Sealing unit for trocars Ø 5mm (20 units)

Single unit price \$20

**Trocar Pins**



**EK051R \$160**  
triangular


**EK045R \$160**  
conical sharp

**EK052R \$160**  
triangular

**EK046R \$160**  
conical sharp

**Trocar Sleeves**

Reusable



**EK010R \$240**

smooth, with stopcock 60 mm




**EK011R \$190**

smooth, without stopcock 60 mm



**EK012R \$290**

threaded, with stopcock 60 mm




**EK013R \$205**

threaded, without stopcock 60 mm




**EK014R \$240**

smooth, with stopcock 110 mm




**EK015R \$190**

smooth, without stopcock 110 mm



**EK016R \$290**

threaded, with stopcock 110 mm



**EK017R \$205**

threaded, without stopcock 110 mm



Ø 10mm

**Sealing Units**

Reusable




**EK083P \$125**  
Sealing unit  
for trocars  
Ø 10/12mm with  
5 mm flap converter



**EK086P \$80**  
Sealing unit  
for trocars  
Ø 10/12mm

**Trocar Pins**



110 mm

**EK062R \$160**  
triangular

**EK056R \$160**  
conical sharp

150 mm

**EK063R \$160**  
triangular

**EK057R \$160**  
conical sharp

**Trocar Sleeves**

Reusable



**EK024R \$230**  
smooth,  
with  
stopcock  
110 mm



**EK025R \$190**  
smooth,  
without  
stopcock  
110 mm



**EK026R \$265**  
threaded,  
with  
stopcock  
110 mm



**EK027R \$205**  
threaded,  
without  
stopcock  
110 mm



**EK028R \$245**  
smooth,  
with  
stopcock  
150 mm





Ø 12mm

**Sealing Units**

Reusable




**EK083P \$125**  
Sealing unit for trocars Ø 10/12mm with 5 mm flap converter



**EK086P \$80**  
Sealing unit for trocars Ø 10/12mm

**Trocar Pins**



**EK072R \$170**  
triangular

**EK066R \$170**  
conical sharp

**EK073R \$170**  
triangular

**EK067R \$170**  
conical sharp

**Trocar Sleeves**

Reusable



**EK034R \$240**  
smooth, with stopcock 110 mm



**EK035R \$190**  
smooth, without stopcock 110 mm



**EK036R \$280**  
threaded, with stopcock 110 mm



**EK037R \$200**  
threaded, without stopcock 110 mm



**EK038R \$270**  
smooth, with stopcock 150 mm

Sealing Units / Consumables

Reusable



**EK083P \$125**  
Sealing unit for trocars Ø 10/12mm with 5 mm flap converter



**EK087P \$370**  
Reducing converter for trocars Ø 10/12mm with 5 mm flap converter (5 units)  
*Single unit price \$80*



**EK086P \$80**  
Sealing unit for trocars Ø 10/12mm



**EK085P \$515**  
Sealing cap for trocars Ø 10/12mm (20 units)  
*Single unit price \$30*




**EK084P \$190**  
Cross slit valve for trocars Ø 10/12mm (20 units)  
*Single unit price \$10*

Trocar for mini-laparotomy (Hasson technique)

Reusable



**EK024R \$225**  
smooth, with stopcock 110 mm



**EK064R \$160**  
blunt (Hasson) 110 mm



**EK098R \$340**  
conical attachment for 10 mm trocar sleeve





Extra Long Cabinet

**\$12,970**

**ARTG  
Listed**

### Andersen Anprolene Ethylene Oxide (EtO) Sterilisation Cabinet

The most gentle, efficient and effective low temperature steriliser available. The Anprolene system offers unique plastic sterilisation bags and individual gas ampoules in two different sizes to match the size of your sterilisation load.

#### Features:

- **SAFE MICRO-DOSE** gas ampoule system – each sterilisation load uses less than 20gr of EtO
- Unmatched sterilisation flexibility to sterilise a wide range of delicate **ENDOSCOPY and VETERINARY** instruments
- Powerful ventilation and purge pumps
- Digital display tracks cycle progress
- Low cost per cycle – between \$25 and \$65 depending upon sterilisation load
- 12 Hour and 24 Hour cycles to choose from
- Extra long cabinet – interior dimensions: 35cm W x 32.5cm H x 74.5cm L / 84L capacity

#### Ideal For:

- Administration Sets
- Arthroscopes & Laparoscopes
- Biopsy & Grasping Forceps
- Bottles
- Catheters
- Cautery Cords & Tips
- Clothes
- Defibrillator Leads & Paddles
- Dermal Punches
- Drill Pieces
- Endotracheal Tubes
- Flexible Endoscopes
- Instruments
- Light Guide Cables
- Gloves
- Needles
- Ophthalmic Instruments
- Oscopes
- Respirators
- Surgical Staplers
- Syringes
- Tubing
- Video Endoscopy Cameras

**Andersen Small and Large Ampoule and Bag Refill Kit**



In the Andersen system, the EtO is supplied in a glass ampoule which is held within a break shield. Gas is delivered to the heart of the load, achieving the highest level of sterilisation for any given amount of gas.

The gas is released by breaking the neck of the ampoule. The EtO boils, giving a very distinctive bubbling and cooling feeling as it does so. This provides very clear, positive feedback that the ampoule has been properly activated.

The packaged items, ready to be sterilised, are placed inside the supplied sterilisation bag. The cabinet holds the space around the sterilisation bag below room pressure. This ensures that any leak is from the room into the cabinet. Any EtO released through the bag is taken to the outside through the exhaust tube at the back of the cabinet.

The bag volumes vary in size, giving great flexibility in the volume, size and shape of items that can be sterilised.



**Small Ampoule and Bag Refill Kit**

**AN7300**

**\$1,290**

consists of:

- 60 x 5cc Ampoules
- 60 x 10 Litre Sterilisation Bags

**Large Ampoule and Bag Refill Kit**

**AN7916**

**\$780**

consists of:

- 14 x 18cc Ampoules
- 14 x 40 Litre Sterilisation Bags
- Probe Strap

**Andersen Dosimeter Chemical Integrator**



Time, temperature and EtO concentration are essential to successful sterilisation. Dosimeter integrators present visual assurance that all these parameters have been fulfilled. Dosimeters change colour in proportion to the dose of sterilant they contact. View results immediately after sterilisation. No laboratory culture is necessary.

Use one Dosimeter per cycle.

**Features:**

- Read results immediately
- Reliable
- Easy to use for every sterilisation load

**AN87**

**\$115**

Box of 25

**Contact Us for Reference Sites**



**ETO Exposure Indicator Tape**

**\$30**  
per roll




**Austvet Endoscopy is committed to:**

- Providing the highest level of service and repair
- Providing the fastest turnaround time
- Maximising our customers' satisfaction and productivity
- Providing training and technical support

**OLYMPUS®**

**AOHUA**



**Don't Take Risks With Your Precious Endoscopy Equipment!**



Our service facility is Olympus and AOHUA approved and all repairs are carried out strictly to manufacturers specification.

**This Guarantees a quality repair**

We are the only Veterinary Supplier in Australia capable of performing the full range of Flexible Endoscopic Repairs In-House.

**This Guarantees a fast turn-around time**



**Don't Take Risks With Your Precious Endoscopy Equipment!**





### Endoscope Cleaning Brushes

Instrument Channel	Working Length	Price
1.0-1.5mm	240cm	<b>\$10</b>
2.0-2.8mm	220cm	<b>\$10</b>
2.0-2.8mm	350cm	<b>\$10</b>

### Matrix Wipes



Unique biofilm removing detergent wipes for pre-cleaning endoscopes and other medical devices. Ready to use and conveniently packaged wipes. Contain MatrixT biofilm removing solution (a patented formulation of surfactants specifically designed for removing biological soils and biofilms). Apertured wipes to ensure sufficient density for removing organic soiling.

**\$10**

### Matrix Biofilm Remover Cleaning Solution



Matrix is a breakthrough in validated cleaning protocols. More effective than enzymatic cleaners, Matrix is the only proven solution designed to remove bacterial biofilms in which both bacteria and other microorganisms may be imbedded on flexible and rigid endoscopes and accessories which have accumulated over time due to insufficient cleaning.

**\$200**  
5 Litre

**Impacts 100% on Viable Bacteria to Provide Absolute Assurance That All Your Endoscopes and Instruments are Optimally Cleaned.**

### Medizyme Multi-Enzyme Cleaning Solution



Medizyme **Multi-Enzyme Cleaner** is totally compatible with all flexible and rigid endoscopes due to its neutral PH. Provides safe and effective cleaning against proteins, carbohydrates, lipids and mucopolysaccharides.

**\$20**  
500ml

Also provides a combination of better detergency and good rinsing ability to ensure the best possible results.

**\$110**  
5 Litre

### OPAL Disinfectant



OPAL™ is an Ortho-Phthalaldehyde (OPA) based **high level liquid disinfectant** that features a reuse period of up to 21 days and a fast turnaround of 6-minute contact time at 20°C. Ready-to-use, no mixing or activation required. Effective against a wide array of micro-organisms. Suitable for use on all flexible and rigid endoscopes and accessories.

**\$150**  
5 Litre

### OPAL Test Strips



Opal® Test Strips are chemical indicators for use in determining whether the concentration of OPA, the active ingredient in OPAL®, is above or below the Minimum Recommended Concentration of 0.35%. It is recommended the OPAL® solution be tested before each disinfection cycle.

**\$60**  
50 Test Strips

- Easy to use - 3 second dip time
- Results in 90 seconds
- Colour reference chart to confirm concentration amount. Large, easy to read indicator pad

### Mr Clear Anti-Fog Solution



Convenient 6ml Anti-Fog solution which eliminates condensation at the distal end of rigid endoscopes/telescopes. Non-abrasive sponge allows for effective cleaning without distal lens damage.

**\$15**



**OLYMPUS®**



**Olympus' CX23 Trinocular microscopes will help you visualise and capture your samples with precise detail, vibrant colours and high resolution.**


The New Olympus CX23 Trinocular microscope offers high optical performance with large field number of 20, a built-in LED light source that provides uniform and stable illumination with low power consumption and an ergonomic design for easy and safe operation.

The CX23 Trinocular microscope can be coupled with a range of cameras that will enable capturing high resolution images faster and easier.

**Specifications:**

- Expected lifetime of LED is over 20,000 hours
- Supplied with 4x, 10x, 40x and 100x Plan Objectives
- Trinocular Head

**\$1,600**



**Introducing  
the NEW  
CX23  
Trinocular  
Microscope**

**MICROSCOPE**

**Accessories:**



Microscope  
Storage/Transport Box

**\$160**



Olympus  
Microscope Dust Cover

**\$12**



Olympus  
Microscope Oil 8ml

**\$15**



**AE**  
AUSTVET ENDOSCOPY

**ANPRO**  
ANDERSEN  
PRODUCTS

ARTG  
Listed



**EtO Sterilisation for the 21<sup>st</sup> Century**

We are also proud to represent Olympus, Aohua, B. Braun & Dr. Fritz

**OLYMPUS**<sup>®</sup>

**AOHUA**

**B | BRAUN**  
SHARING EXPERTISE

**Dr. Fritz**<sup>®</sup>  
GmbH