

Standardization
with Specialization



Joint Replacements

Trauma

Spine

Micro Implants

Orthobiologics

Instruments

Interventional Pain

Navigation

Endoscopy

Communications

Patient Handling Equipment

EMS Equipment

The information presented in this brochure is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Products referenced with ™ designation are trademarks of Stryker.
Products referenced with ® designation are registered trademarks of Stryker.

Literature Number: 9100-927-002 Rev. None
UnDe/NAC

Copyright © 2004 Stryker
Printed in USA

Stryker Instruments
4100 East Milham Avenue
Kalamazoo, MI 49001 USA
t: 269 323 7700 f: 800 999 3811
toll free: 800 253 3210

www.stryker.com

High Speed Drill Systems



Specialty-driven handpieces

Increased power and versatility

Complete OR standardization

Standardization

Solutions for the OR

Simplify the Operating Room

- One line of products for all needs
- Reduce bur and blade inventory
- Streamline staff education

Quality Stryker Service

- Free loaner program
- Original Stryker parts and labor
- Preventative maintenance programs
- One source, one solution

Financial Options

- Stabilize operational costs
- Total financing
- Customized solutions

pneumatic electric

Maestro Drill

Spine/Neuro/ENT

5400-200

- More power at lower PSI
- Double swivel mechanism
- Elite telescoping burs
- Decreased noise and increased reliability

Elite Burs

- One length of bur will fit every PD/SD and Elite attachment to simplify and reduce inventory
- Numbered safety marks help surgeon and staff adjust bur exposure

Saber Drill

ENT/Spine

5400-120

- Elite telescoping burs
- PD/SD Series attachments

Maestro Footpedal

5400-208

- Variable speed control
- Small footprint



Diffuser

5400-006 (10/pkg)

Oil Cartridge

5400-005 (10/pkg)

- Separate oil cartridge and diffuser for more consistent lubrication

U Drill

Spine/Neuro

5400-100

- Electronic torque feedback
- Tapered attachments

UHT Drill

Neuro

5400-110

- Fixed and steering duraguards
- Elite attachments



CORE Console

5400-050

- Operates two handpieces simultaneously
- Integrated irrigation pump
- Intuitive color touch screen display

stryker solutions